



Agilent NMR Systems AutoTest

User Guide



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2013

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G7446-33301

Edition

Rev. A, March 2013

Printed in USA

Agilent Technologies, Inc.
5301 Stevens Creek Boulevard
Santa Clara, CA 95051 USA

Warranty

The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or sub-contract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14

(June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

Contents

1 Introduction

About This Guide 10

2 AutoTest Interface

Introduction 12

AutoTest Configuration Tab 14

Selecting Test Packages 16

AutoTest Test Library Tab 16

General 16

2H tests 17

1H tests 17

13C tests (using channel 2 hardware - normal cabling) 18

13C tests (using channel 3 hardware - must re-cable) 18

Shaped pulses 18

Standard gradient tests (Z) 19

Standard gradient tests (X) 19

Standard gradient tests (Y) 19

Power-limited 13C pulse tests (channel 2) 19

Power-limited 13C pulse tests (channel 3) 20

Power-limited 13C/N15 decoupling efficiency tests 20

Power-limited 13C/N15 decoupling noise and stability tests 20

13C N15 decoupling installation tests for cryogenic probe 20

Extended gradient tests (Z) 21

Extended gradient tests (X) 21

Extended gradient tests (Y) 21

Probe overnight tests 22

N15 tests (using channel 3 hardware - cable to N15 probe port) 22

N15 tests (using channel 2 hardware - cable to N15 probe port) 22

AutoTest sample with 100 mM NaCl 22

AutoTest sample with 250 mM NaCl 23

Power-limited 13C Pulse / decoupling tests (Channel 2: approx 10
microsecond 13Cpw90) 23

Probe overnight tests (approx 10 microsecond 13Cpw90) 23

History Tab 24

Test Report 25

AutoTest Settings 26

3 AutoTest Operation

Getting Started	34
AutoTest sample	34
AutoTest for ATP (acceptance test procedure) and running it the first time	34
System setup	34
AutoTest Setup	36
Start AutoTest	36
Enter system information	36
Configure AutoTest	38
Output configuration	38
Enable tests	39
Select test packages	40
Run AutoTest	42
Standard Tests Performed by AutoTest	42
Automated console tests	42
Automated tests with shaped RF	43
Automated decoupling performance tests	43
Automated 90 degree pulse width calibrations (PW90)	44
RF homogeneity tests	44
Gradient calibrations and performance tests	44

4 AutoTest Administration

Saving Data and FID Files from Previous Runs	48
Creating Probe-Specific Files	48
AutoTest Directory Structure	49
data directory	50
data.failed directory	50
parameters directory	51
reports directory	51
texts directory	51
history directory	51
atdb directory	52
AutoTest Macros	54
ATglobal macro	54
ATstart macro	57
ATnext macro	57
ATxxx macros	57
ATcxxx macros	58
ATrecord macro	58

5 AutoTest - Test Reference

RF Performance Test Descriptions	62
Pulse Stability and Sensitivity	62
Cancellation test	63
Phase stability (13 deg. phase error) test	63
Pulse turn-on time	64
Attenuator linearity	64
Attenuator linearity at reduced power	64
Linear modulator linearity tests	65
Linear modulator linearity tests with attenuators set to low power	65
Pulse shape test-DANTE	65
Phase switch settling time	66
RF homogeneity	66
Receiver test	67
Image rejection test	68
Quadrature phase selection	68
Shaped Pulse Test Descriptions	69
Gaussian-shaped pulse excitation	69
Gaussian 90 degree pulse stability	69
Gaussian 13 degree phase error	70
Gaussian SLP 13 degree phase error (phase-ramped gaussian pulses)	70
Shaped pulse settability	70
Shaped pulse test - rectangular, gaussian, and EBURP-1	70
Shaped pulse test-constant bandwidth for a variety of shapes	71
Shaped pulse scalability	71
13C Test Descriptions	72
13C 90 degree pulse width calibration	72
X-Phase modulation decoupling profiles	72
13C-Adiabatic decoupling profiles	72
Gradient Tests Descriptions	74
Gradient profile	74
Field recovery stability	74
Field recovery	74
Gradient echo stability	75
Gradient effect on cancellation test	75
CPMG T_2	75
13C Power-Limited Pulse Tests	77
13C pw90 (15 microsecond) and lowband amplifier compression	77
13C pw90 (10 microsecond) and lowband amplifier compression	77

13C Power for pw90=15.0 microsecond and lowband amplifier compression	77
13C Power for pw90=10.0 microsecond and lowband amplifier compression	78
13C RF homogeneity at limited power (15 microsecond)	78
13C RF homogeneity at limited power (10 microsecond)	78
13C RF homogeneity at limited power (pw90=15.0 microsecond)	78
13C RF homogeneity at limited power (pw90=10.0 microsecond)	79
13C RF amplifier linearity at limited power (pw90=15 microsecond)	79
13C RF amplifier linearity at limited power (pw90=10 microsecond)	79
13C/15N power-limited decoupling tests	81
13C/15N decoupling noise (FID) at user-selected decoupling power	81
13C/15N decoupling noise (FID) as function of decoupling power	81
Sensitivity as function of 13C decoupling power	82
Sensitivity as function of 13C decoupling power in presence of 15N decoupling	82
Sensitivity as function of 15N decoupling power	83
Sensitivity as function of 15N decoupling power in presence of 13C decoupling	84
13C decoupling using phase modulation	84
13C decoupling using adiabatic decoupling	85
Decoupler heating	85
Amplitude stability in the presence of 13C decoupling	85
Amplitude stability in the presence of 15N decoupling	85
Amplitude stability in the presence of combined 15N and 13C decoupling	85
Decoupled methanol amplitude stability using a 6 kHz 13C decoupling RF field	86
Decoupled methanol amplitude stability using a 3kHz 13C decoupling RF field	86
Comparison of single-pulse versus HSQC decoupled methanol	86
Comparison of single-pulse versus HSQC decoupled methanol (13C pw90=10 microsecond)	87
15N Power-Limited Pulse and Decoupling Tests	88
15N pw90 and lowband amplifier compression	88
15N RF amplifier linearity at limited power	88
Acetonitrile 15N-HMQC amplitude stability	88
Acetonitrile 15N-HMQC decoupler heating	88
2H Pulse and Decoupling Tests	89
2H pw90 using channel 2	89
2H pw90 using channel 3	89
2H pw90 using channel 4	90
2H pw90 using lock/decoupler	90

^2H 90 degree pulse amplitude stability and sensitivity	90
^2H 90 degree pulse amplitude stability using channel 3	90
^2H 90 degree pulse amplitude stability using channel 4	90
^2H 90 degree pulse amplitude stability using lock/decoupler	92
^2H spinlock test	92
Installation Tests for Cryogenic Probes	92
15-minute conditioning (ATcryo15)	93
30-minute conditioning (ATcryo30)	93
60-minute conditioning (ATcryo60)	93
90-minute conditioning (ATcryo90)	93
Other Test Descriptions	93
Heating during spin lock test	93
Lock tests	94
Spectral purity test	94
Variable temperature test	94
Small-angle phase shift test	95
Tests Using Salty Sample	95
^1H RF homogeneity	95
Pulse stability	95
Heating during spin lock test	96
Decoupler heating	96



1

Introduction

About This Guide 10

Agilent MR systems provide unmatched productivity for diverse chemical applications by combining easy-to-use software with outstanding performance. Push-button experiments, along with straightforward processing and data export capabilities, make the MR the best choice for compound detection, quantification and structure confirmation.



About This Guide

The following chapters are included in this guide:

[Chapter 2](#), “AutoTest Interface”

[Chapter 3](#), “AutoTest Operation”

[Chapter 4](#), “AutoTest Administration”

[Chapter 5](#), “AutoTest - Test Reference”



2 AutoTest Interface

Introduction	12
AutoTest Configuration Tab	14
Selecting Test Packages	16
AutoTest Test Library Tab	16
History Tab	24
Test Report	25
AutoTest Settings	26



Introduction

The AutoTest interface is comprised of the **AutoTest** Settings window and the **AutoTest Run** window. Both windows are accessed from the **Tools > Standard Calibrations Experiments** menu.

To open the **AutoTest Run** window, select **Tools > Standard Calibration Experiments > Start Autotest**.

The AutoTest window is a graphical user interface for configuring and running tests. It features a tabbed interface with the following sections:

- Configuration:** Contains input fields for:
 - Operator Name: vnmr1
 - Console type: VNMR5
 - Console S/N: 3415
 - Data Directory: /home/vnmr1/vnm
 - Probe: HCN
 - Coil size(mm): 16
 - Max 13C Pulse Power: 60
 - Max 13C Dec.Power: 46
 - Max 15N Pulse Power: 60
 - Max 15N Dec.Power: 44
 - 1H Pulse Power: 62
 - Nominal 1H pw90: 7
- Test Library:** Contains a list of test categories with checkboxes:
 - ☒ Auto Plotting
 - ☐ Print Parameters
 - ☐ Graph History
 - ☒ Process During Acq.
 - ☒ Plot to file (PostScript printer only)
 - ☒ Repeat Until Stopped
 - ☒ Enable VT tests
 - ☒ Enable C13 tests
 - ☒ Enable Lock tests
 - ☒ Enable Gradient tests
 - ☐ All Standard Tests
 - ☐ Signal Detection
 - ☐ All Channel 1 (1H) Tests
 - ☐ Channel 1 Non-Shaped Tests
 - ☐ Channel 1 Shaped Tests
 - ☐ Z-Gradient Tests
 - ☐ X-Gradient Tests
 - ☐ Y-Gradient Tests
 - ☐ All Standard Tests (Power Limited)
 - ☐ 13C Tests (channel 2)
 - ☐ 13C Tests (channel 3 -recable)
 - ☐ 13C Tests using Above Power Limits (Channel 2)
- History:** (Empty)
- Test Report:** (Empty)

At the bottom of the window are four buttons: **Begin Test**, **Clear Selection**, **Print Graph**, and **Exit**.

Figure 1 The AutoTest window

To open the **AutoTest Settings** window, select
Tools > Standard Calibration Experiments > Autotest settings.




Figure 2 The AutoTest setting window

AutoTest Configuration Tab

The **Configuration** tab on the **AutoTest Run** window is displayed when **AutoTest** is started. All the information required to set up AutoTest, identify the system and probe, power levels, and other information are entered in the fields shown on the **Configuration** tab.

To open **AutoTest Run** window and view the **Configuration** tab select **Tools > Standard Calibration Experiments > Start Autotest**.

The screenshot shows the 'Auto Test' window with the 'Configuration' tab selected. The window is divided into several sections:

- System information fields:** This section contains input fields for:
 - Operator Name: vnmr1
 - Console type: VNMRS
 - Console S/N: 3415
 - Data Directory: /home/vnmr1/vnm
 - Probe: HCN
 - Coil size(mm): 16
 - Max 13C Pulse Power: 60
 - Max 13C Dec.Power: 46
 - Max 15N Pulse Power: 60
 - Max 15N Dec.Power: 44
 - 1H Pulse Power: 62
 - Nominal 1H pw90: 7
- Configure AutoTest:** This section contains checkboxes for:
 - Auto Plotting
 - Print Parameters
 - Graph History
 - Process During Acq.
 - Plot to file (PostScript printer only)
 - Repeat Until Stopped
 - Enable VT tests
 - Enable C13 tests
 - Enable Lock tests
 - Enable Gradient tests
- Click Arrows to Navigate between Test Pages:** This section contains two arrows pointing left and right.
- Available AutoTests (page 1 of 2):** This section contains a list of test categories:
 - All Standard Tests
 - Signal Detection
 - All Channel 1 (1H) Tests
 - Channel 1 Non-Shaped Tests
 - Channel 1 Shaped Tests
 - Z-Gradient Tests
 - X-Gradient Tests
 - Y-Gradient Tests
 - All Standard Tests (Power Limited)
 - 13C Tests (channel 2)
 - 13C Tests (channel 3 -recable)
 - 13C Tests using Above Power Limits (Channel 2)

At the bottom of the window are four buttons: **Begin Test**, **Clear Selection**, **Print Graph**, and **Exit**.

Figure 3 Configuration tab - Tests page 1 of 2

Auto Test

Configuration Test Library History Test Report

Operator Name: vnmr1 Max 13C Pulse Power: 60

Console type: vnmrs Max 13C Dec.Power: 45

Console S/N: 41005 Max 15N Pulse Power: 50

Data Directory: /home/vnmr1/vnm Max 15N Dec.Power: 34

Probe: HCN 1H Pulse Power: 54

Coil size(mm): 18 Nominal 1H pw90: 6

☒ Auto Plotting
 ☐ Enable VT tests

☐ Print Parameters
 ☒ Enable C13 tests

☐ Graph History
 ☐ Enable Lock tests

☒ Process During Acq.
 ☒ Enable Gradient tests

☒ Repeat Until Stopped

☐ 13C Tests using Above Power Limits (Channel 3 -recable)

☐ 1H/13C Tests using 100mM NaCl Sample

☐ 1H/13C Tests using 250mM NaCl Sample

Figure 4 Configuration tab - Tests page 2 of 2

Selecting Test Packages

The **Configuration** tab has two pages of test packages. These pages provide access to a standard test package by clicking on **All Standard Tests** or individual test packages. The individual tests used in each test package are accessed by clicking on the **Test Library** tab. Refer to “[AutoTest Test Library Tab](#)” on page 16 for the list packages and the tests in each package.

AutoTest Test Library Tab

AutoTest has a comprehensive set of experiments that can be run alone or several experiments can be run as a package. Individual tests are accessed from the **Test Library** tab. Some individual tests are used in multiple packages and therefore appear several times.

- Click the right arrow button to move to the next page of tests, or the left arrow to move back one page. The right or left arrow is no longer visible when the first or last page is displayed.
- Mark the check box next to an experiment to enable the experiment.

To display a list of specific tests select the **Test Library** tab at the top of the **AutoTest** window. Use the right and left arrow buttons to view groups of individual tests organized by hardware; for example, channel 1 tests, z-axis gradient tests, ¹³C tests, and so on.

Select one or more of the check boxes in any of these to display the tests to be performed, in the order selected. To clear a test, deselect the check box. Any of these tests may be grouped with the group test(s) specified in the **Configuration** panel.

To initiate an experiment, select the **Configuration** tab and click the **Begin Test** button. Experiments in the test library are grouped under the following headings.

General

- H1 RF Homogeneity
- High-Frequency Amplifier Compression
- T1 Determination

- Receiver Gain
- DSP Sensitivity Improvement
- Quadrature Image
- Spectral Purity Test
- Phase-Cycle Cancellation
- Phase-Cycle Cancellation vs. Recycle Time
- Spinlock Heating Test
- variable Temperature Test
- Lock Power and Gain Text
- Gradient Mapping and Shimming
- RF Field Mapping

2H tests

- 2H PW90 on Lock Coil (tn = 1k)
- 2H Stability / Sensitivity (tn = 1k)
- 2H spinlock Test (tn = 1k)
- 2H PW90 on Lock Coil (Channel 3, Recable) For Inova only
- 2H Stability (Channel 3, Recable) For Inova only
- 2H PW90 on Lock Coil (Channel 4, Recable) For Inova only
- 2H Stability (Channel 4, Recable) For Inova only
- 2H PW90 on Lock Coil (Lock / Decoupler, Recable) For Inova only
- 2H Stability (Lock/Decoupler, Recable) For Inova only
- 2H PW90 on Lock Coil (Channel 5, Lock/Decoupler, Recable) For Inova only
- 2H Stability (Channel 5, Lock/Decoupler, Recable) For Inova only

1H tests

- H1 pw90
- 90/10 Degree, 1 μ s Pulse Stability/Sensitivity
- 30 Degree Pulse Stability
- Phase Stability
- Quadrature Phase Shift

- Small-Angle Phase Shift
- Phase Switch Timing
- Pulse Turn-on Time
- Dante Pulse Train
- Attenuator Test
- Linear Modulator Test

13C tests (using channel 2 hardware - normal cabling)

- 13C PW90 and Low-Band Amplifier Compression
- 13C RF Homogeneity
- 13C Phase Modulation Decoupling Profiles
- 13C Adiabatic Decoupling Profiles
- Methanol Amplitude Stability using 3kHz 13C Decoupling
- Methanol Amplitude Stability using 6kHz 13C Decoupling
- 13C Decoupling Heating Test

13C tests (using channel 3 hardware - must re-cable)

- 13C PW90 and Low-Band Amplifier Compression
- 13C RF Homogeneity
- 13C Phase Modulation Decoupling Profiles
- 13C Adiabatic Decoupling Profiles
- 13C Decoupling Heating Test

Shaped pulses

- Gaussian Pulse Excitation
- Gaussian 90 Degree Stability
- Gaussian Phase Stability
- GaussianSLP Phase Stability
- Gaussian Power / Pulse Width Array
- Rectangular vs. Gauss vs. EBURP1
- Pbox Shapes
- Shaped Pulse Amplitude Scaling

Standard gradient tests (Z)

- Z-Gradient Calibration
- Z-Gradient (Field) Recovery Stability
- Z-Gradient (Field) Recovery Rect (30G/cm)
- Z-Gradient (Field) Recovery Sine (30G/cm)
- Z-Gradient (Field) Recovery Rect (10G/cm)
- Z-Gradient (Field) Recovery Sine (10G/cm)
- Phase-cycle cancellation after Gradient
- Z-Gradient Echo Stability (30G/cm)
- Z-Gradient Echo Stability (10G/cm)
- Z-Gradient CPMGT2

Standard gradient tests (X)

- X-Gradient Calibration
- X-Gradient (Field) Recovery Stability
- X-Gradient (Field) Recovery Rect (10G/cm)
- X-Gradient (Field) Recovery Sine (10G/cm)
- X-Gradient Echo Stability (10G/cm)
- X-Gradient CPMGT2

Standard gradient tests (Y)

- Y-Gradient Calibration
- Y-Gradient (Field) Recovery Stability
- Y-Gradient (Field) Recovery Rect (10G/cm)
- Y-Gradient (Field) Recovery Sine (10G/cm)
- Y-Gradient Echo Stability (10G/cm)
- Y-Gradient CPMGT2

Power-limited ¹³C pulse tests (channel 2)

- ¹³C PW90 and LB Amplifier Compression
- ¹³C RF Homogeneity
- ¹³C PW90 / Fine Power for 15.0 μ s 90
- ¹³C RF Homogeneity for 15.0 μ s 90
- ¹³C pw90 vs. power: Amplifier Linearity

Power-limited ¹³C pulse tests (channel 3)

- ¹³C PW90 and LB Amplifier Compression
- ¹³C RF Homogeneity
- ¹³C pw90 vs. power: Amplifier Linearity

Power-limited ¹³C/¹⁵N decoupling efficiency tests

- ¹³C Phase Modulation Decoupling Profiles
- ¹³C Adiabatic Decoupling Profiles
- ¹³C WURST2 / WRUST40 Decoupling Profiles
- ¹³C Adiabatic, Waltz, No Decoupling 1D Comparisons
- ¹³C Decoupling 1D Comparisons
- ¹³C 14kHz Adiabatic Decoupling Intensity / Linewidth Stability
- ¹³C 6kHz Adiabatic Decoupling Intensity / Linewidth Stability
- ¹³C 140ppm Adiabatic Decoupling Intensity / Linewidth Stability
- ¹³C Decoupling Heating Test
- ¹³C Decoupling Line Broadening Test

Power-limited ¹³C/¹⁵N decoupling noise and stability tests

- ¹³C / ¹⁵N Decoupling Noise in FID at Max Power
- ¹³C / ¹⁵N Decoupling Noise in FID vs. Power
- Methanol Amplitude Stability using 3kHz ¹³C Decoupling
- Methanol Amplitude Stability using 6kHz ¹³C Decoupling
- H₂O Amplitude Stability with ¹³C Decoupling
- H₂O Amplitude Stability with ¹⁵N Decoupling
- H₂O Amplitude Stability with ¹³C and ¹⁵N Decoupling
- H₂O SN vs. ¹³C Decoupling Power
- H₂O SN vs. ¹³C Decoupling Power (w/¹⁵N dec)
- H₂O SN vs. ¹⁵N Decoupling Power
- H₂O SN vs. ¹⁵N Decoupling Power (w/¹³C dec)

¹³C ¹⁵N decoupling installation tests for cryogenic probe

- 15 Minute Conditioning Alternated with CryoNoise Test

- 30 Minute Conditioning Alternated with CryoNoise Test
- 90 Minute Conditioning Alternated with CryoNoise Test

Extended gradient tests (Z)

- Z-Gradient Calibration
- Z-Gradient (Field) Recovery Stability
- Z-Gradient (Field) Recovery Rect (30G/cm)
- Z-Gradient (Field) Recovery Rect (30G/cm, lb=1)
- Z-Gradient (Field) Recovery Sine (30G/cm)
- Z-Gradient (Field) Recovery Rect (10G/cm)
- Z-Gradient (Field) Recovery Sine (10G/cm)
- Phase-cycle cancellation after Gradient
- Z-Gradient Echo Stability (30G/cm)
- Z-Gradient Echo Stability (10G/cm)
- Z-Gradient CPMGT2
- Z-Gradient CPMGT2 (30G/cm)

Extended gradient tests (X)

- X-Gradient Calibration
- X-Gradient (Field) Recovery Stability
- X-Gradient (Field) Recovery Rect (10G/cm)
- X-Gradient (Field) Recovery Rect (maxG/cm)
- X-Gradient (Field) Recovery Sine (10G/cm)
- X-Gradient Echo Stability (10G/cm)
- X-Gradient CPMGT2
- X-Gradient CPMGT2 (20G/cm)

Extended gradient tests (Y)

- Y-Gradient Calibration
- Y-Gradient (Field) Recovery Stability
- Y-Gradient (Field) Recovery Rect (10G/cm)
- Y-Gradient (Field) Recovery Rect (maxG/cm)
- Y-Gradient (Field) Recovery Sine (10G/cm)
- Y-Gradient Echo Stability (10G/cm)
- Y-Gradient CPMGT2

- Y-Gradient CPMGT2 (20G/cm)

Probe overnight tests

- ^{13}C HSQC Stability
- ^{13}C HSQC Stability (power limited)
- ^{13}C , ^{15}N HSQC-NOESY (power limited)
- ^{13}C HSQC 1D vs. s2pul (power limited)
- Methanol HSQC Stability vs. ^{13}C Pulse Power (coupled)
- Methanol HSQC Stability vs. ^{13}C Pulse Power (decoupled)
- Methanol HSQC Stability vs. ^{13}C Decoupling Power (decoupled)

N15 tests (using channel 3 hardware - cable to N15 probe port)

- N15 PW90 and Low-Band Amplifier Compression
- N15 HMQC Stability
- N15 pw90 vs. Power: Amplifier Linearity
- N15 Decoupling Heating Test
- N15 Decoupling Heating Test with Power Limit

N15 tests (using channel 2 hardware - cable to N15 probe port)

- N15 PW90 and Low-Band Amplifier Compression
- N15 HMQC Stability
- N15 pw90 vs. Power: Amplifier Linearity
- N15 Decoupling Heating Test
- N15 Decoupling Heating Test with Power Limit

AutoTest sample with 100 mM NaCl

- H1 RF Homogeneity (with 100mM NaCl)
- H1 pw90 (with 100mM NaCl)
- 90- and 10-Degree Stability and sensitivity (with 100 mM NaCl)
- Spinlock Heating Test
- ^{13}C Decoupling Heating Test

AutoTest sample with 250 mM NaCl

- H1 RF Homogeneity (with 250mM NaCl)
- H1 pw90 (with 100mN NaCl)
- 90- and 10-Degree Stability and Sensitivity (with 250 mM NaCl)
- Spinlock Heating Test
- ^{13}C Decoupling Heating Test

Power-limited ^{13}C Pulse / decoupling tests (Channel 2: approx 10 microsecond $^{13}\text{Cpw90}$)

- ^{13}C PW90 and LB Amplifier Compression
- ^{13}C RF Homogeneity
- ^{13}C PW90 for 10 μs 90
- ^{13}C RF Homogeneity for 10 μs 90
- ^{13}C pw90 vs. Power: Amplitude Linearity
- Methanol Amplitude Stability using 3kHz ^{13}C Decoupling
- Methanol Amplitude Stability using 6kHz ^{13}C Decoupling
- ^{13}C Decoupling Heating Test

Probe overnight tests (approx 10 microsecond $^{13}\text{Cpw90}$)

- ^{13}C HSQC Stability ($\sim 10 \mu\text{s}$ $^{13}\text{Cpw90}$)
- ^{13}C HSQC Stability ($\sim 10 \mu\text{s}$ $^{13}\text{Cpw90}$ Power Limited)
- ^{13}C , ^{15}N HSQC-NOESY ($\sim 10 \mu\text{s}$ $^{13}\text{Cpw90}$ Power Limited)
- ^{13}C HSQC 1D vs. s2pul ($\sim 10 \mu\text{s}$ $^{13}\text{Cpw90}$ Power Limited)
- Methanol HSQC Stability versus ^{13}C Pulse Power (Coupled)
- Methanol HSQC Stability versus ^{13}C Pulse Power (Decoupled)
- Methanol HSQC Stability versus ^{13}C Decoupling Power (Decoupled)

History Tab

Selecting the **History** tab at the top of the **AutoTest** window provides a graphical output of the history files accumulated after several AutoTest runs.

To select a history file, scroll to the name of the file. Use the left arrow and right arrow to rapidly step through the list.




Figure 5 AutoTest - History tab

Buttons are provided to specify a graphical or text output of the history file. If more than one result is contained in the history file, a menu of choices is displayed to the right of the display of results.

Values of results that exceed specified limits are displayed in red in the graphical output and in color in the text output. If a history file has multiple results with specifications, all displayed points for that run will indicate failure if only one has failed. Thus, a particular display may show a failure

even though the result is within specification. If this is true, select the other entries in the list to show which result had actually failed.

Printed graphs may be obtained by selecting the button at the bottom of the screen. A full set of small graphs is plotted automatically after an AutoTest run when the **Graph History** check box is selected in the **Configuration** display (if automatic plotting is requested).

Test Report

Selecting the **Test Report** tab at the top of the **AutoTest** window allows viewing the `~/vnmrsys/autotest/atrecord_report` file, a report on the results for the current AutoTest run. The `~/vnmrsys/autotest/atrecord_report` file is one of two files contains test results (both files are stored in the `autotestdir` directory, typically `~/vnmrsys/autotest/`):

- The `REPORT` file is more compact and is automatically printed at the end of the run.
- The `atrecord_report` has the same results, but in a format that is similar to the history file format. In addition, it indicates a Fail status if any of the results in a history file line is out of bounds relative to the upper and lower limits stored in the `~/vnmrsys/autotest/atdb/at_specs_table` file.

Scroll the `atrecord_report` report to view all of the results. All failures are indicated. Tests that have passed defined upper and lower limits and tests for which there are no defined upper and lower limits have no Pass/Fail status indicated. Select the **Fail Only** option to see a rapid indication of any failures. Select the **Print Report** button, at the bottom of the display, to print the report.




Figure 6 Test Report tab of the AutoTest Window

Failed experimental data is stored, with a date-stamp, in a `data.failed` directory, which is created in the `~/vnmrsys/autotest` directory (if needed). If tests are repeated, new date-stamped files are created at each failure, permitting later examination of the failed experiments. In addition, a file `FAILREPORT` is stored in the same directory and it is printed at the same time as the `REPORT` file. This report shows the history file entry for all failures and the time of the failure. This report is also stored with date-stamp in the `autotestdir/reports` directory when a new **AutoTest** run is begun.

AutoTest Settings

Access the **AutoTest Settings** window by selecting **Tools>Utilities > Standard Calibration Experiments > AutoTest Settings**. **AutoTest** updates the values for each of these settings for which new calibrations are specified when **AutoTest** is run. These results are stored in `~/vnmrsys/global` as they are determined (see “[ATglobal macro](#)” on page 54). There are seven tabs in the **AutoTest Settings** window.




Figure 7 Autotest Settings - Genera




Figure 8 Autotest Settings - 13C Channel 2




Figure 9 Autotest Settings - N15 Channel 3/Channel 2




Figure 10 Autotest Settings - 13C Channel 3




Figure 11 Autotest Settings - 800 MHz




Figure 12 Autotest Settings - H2 Obs/Decouple




Figure 13 Autotest Settings - System



3 AutoTest Operation

Getting Started	34
AutoTest Setup	36
Run AutoTest	42
Standard Tests Performed by AutoTest	42



Getting Started

AutoTest sample

AutoTest uses 0.1% ^{13}C -enriched methanol in 1% $\text{H}_2\text{O}/99\%$ D_2O . More recent samples have 0.1% ^{15}N -enriched acetonitrile added, as well. The sample is doped with gadolinium chloride at a concentration of 0.30 mg/ml, which produces a ^1H T_1 relaxation time of about 50 to 75 ms. The resulting line width is considerably larger than the magnet-determined line width because of the paramagnetic relaxation contribution, minimizing any dependence on shimming skill and permitting rapid collection of data.

AutoTest for ATP (acceptance test procedure) and running it the first time

The ATP is run using a room temperature, 5-mm, Z-axis, Indirect Detection probe (not a Cold Probe).

The first time AutoTest is run for a specific probe, and for the ATP, the full set of ATP tests must be run.

Do the following:

- 1 Follow the procedures in “[AutoTest Setup](#)” on page 36.
- 2 Enable the typical and ATP outputs listed in “[Output configuration](#)” on page 38.
- 3 Enable the typical and ATP tests listed in “[Enable tests](#)” on page 39.
- 4 Select the test package ALL, see “[Select test packages](#)” on page 40.
- 5 Make a backup copy of the final ATP results.

System setup

- 1 Use the sample described in “[AutoTest sample](#)” on page 34.
- 2 Set up all elements of the rf system (transmitters, linear modulators, rf attenuators, amplifiers, receivers, and probes) in their standard configuration.

- 3 If the system has a PFG (pulsed field gradients) accessory installed, do the following
 - a Turn on the gradient amplifier.
 - b Set pfgon correctly for the number of gradients available as follows:
 - Click the **Tools>Utilities** menu.
 - Select **System Settings**.
 - Set the active gradients.
 - c Allow sufficient time for stabilization before running AutoTest. AutoTest will calibrate the gradients.
- 4 Set the VT to 25° C.

Allow the temperature of the sample to regulate and equilibrate. Specifications are determined at 25° C; normal day-to-day AutoTest runs may be done at any temperature, but some results (such as T₁ or T₂ are temperature-dependent).
- 5 Tune the probe, and lock on the D₂O resonance.
- 6 Shim the field on the sample to give a non-spinning half-height that is dominated by the paramagnetic relaxation agent. This should require little time since the H₂O line is quite broad. The methanol doublet is sharper and can be used for a more sensitive test of lineshape and resolution. The acetonitrile will show a doublet in a very well-shimmed system. Resolution-enhancement processing should show a doublet.
- 7 Use a tpwr value necessary to produce a ¹H 90° pulse width that is within (8 to 10 μs) and does not cause probe arcing.

AutoTest uses this tpwr value, determines the ¹H 90° pulse width, and calculates the amplifier compression at this tpwr value.
- 8 If you want to save the results of a previous AutoTest run, rename the history and data directories.
- 9 Select the proper printer from the menu (**File > Printers...**) if hardcopy or Postscript plots are desired.

AutoTest Setup

Set up AutoTest in the following order:

- 1 “Start AutoTest”
- 2 “Enter system information”
- 3 “Configure AutoTest”

CAUTION

Remove any files used for previous versions of AutoTest from your `vnmr`sys directories, particularly any `~/vnmr`sys/seqlib/AT* or `~/vnmr`sys/maclib/AT* files. Remove AT*.DEC, gauss32.RF, gauss.RF and eburp1.RF from `~/vnmr`sys/shapelib if present. Make sure that `~/vnmr`sys/maclib/autotest is not present.

Start AutoTest

- 1 Select **Tools>Utilities > Standard Calibration Experiments > Start AutoTest** from the menu bar.
- 2 The AutoTest widow is displayed.

The system first checks for an autotest directory in the user's current `vnmr`sys directory. If no autotest directory is present, `~/vnmr`sys/autotest is automatically created and subdirectories are copied from `/vnmr/autotest`. These directories include parameter, data, history, and database atdb directories.

Enter system information

At the top of the Configuration display are fields for entering user, console, and probe identification as well as the nominal `tpwr` and `pw90` values for ^1H .

The `tpwr` and `pw90` values are used as starting points and should be entered the first time AutoTest is run. The ^1H `pw90` value should be correct to within $\pm 30\%$. The power level entered is used for most of the tests; so it should be a level that is used in normal day-to-day work. For most probes this would be a power level sufficient to get a 10 μs `pw90`.

Running at the upper range of `tpwr` may result in a situation where the amplifier is in compression. This means that although the `pw90` may be shorter at `tpwr` = 63 than `tpwr` = 57, the `pw90` at the former power may not be a factor of two shorter. If a 100-watt high-frequency amplifier is

present, the t_{pwr} value should be set so that 1H pw90 is 10 μs at that power. In any case, the amplifier compression at the entered t_{pwr} is determined, with the assumption that $t_{pwr} - 12$ dB is a setting where there is linear behavior.

The length (in mm) of the active window in the receiver coil is also specified. This is normally 16 mm, unless the probe has a hexagonal base (Unibody-style), in which case the coil may have an 18 mm length (normally true for indirect and triple-resonance liquids probes).

Some probes (very high field and cryogenic) require more careful power control. Power control for ^{13}C and ^{15}N pulses and decoupling is available within this panel. These values are used to set upper limits on power (attenuator) settings for ^{13}C and ^{15}N (using channels 2 and 3) pulses and decoupling for those tests that have power limits. These are primarily the decoupling noise tests. The ^{13}C pw90 calibration is normally done for a power level giving a nominal 15 μs pulse and this power level is found by experiment. This power value is limited by the above values if the tests are labeled “power limited”. If the probe used has no ^{15}N port, set the maximum pulse power to zero. The relevant macros check for a zero value and skip ^{15}N pulses or decoupling for this case. AutoTest macros beginning with ATC have power limits internally coded.

Enter the following information in their respective data fields:

Table 1 System information

Field label	Field description
Operator name:	Accept the default, the log in name of the current LINUX/UNIX user that has started VnmrJ 3, or type the desired operator name.
Console type:	Enter a description of the console type, such as, 3chan500.
Console S/N:	Enter the console serial number.
Data directory:	Accept the default data path or provide an alternate path, beginning with the root (/).
Probe:	Enter a description of the probe or use the probe's serial number.
Coil size (mm):	Enter the coil size if the probe is equipped with a gradient or gradients coil.

Table 1 System information (continued)

Field label	Field description
Max ^{13}C Pulse Power:	Enter the maximum power (attenuator settings) for any ^{13}C pulse.
Max ^{13}C Dec. Power:	Enter the maximum power (attenuator setting) for any ^{13}C decoupling.
Max ^{15}N Pulse Power:	Enter the maximum power (attenuator settings) for any ^{15}N pulse.
Max ^{15}N Dec. Power:	Enter the maximum power (attenuator setting) for any ^{15}N decoupling.
^1H Pulse Power:	Enter a t_{pwr} for a nominal ^1H pw90 of 10 μs .
Nominal ^1H pw90:	Enter a nominal ^1H pw90 at the above power.

Configure AutoTest

Below the system information entry fields are check boxes for configuring the output and enabling various tests. Complete the AutoTest Configuration in the following order. See:

- 1 “Output configuration” on page 38
- 2 “Enable tests” on page 39
- 3 “Select test packages” on page 40

Output configuration

Check boxes for configuring the output are:

Table 2 Output configuration

Option	Typical and ATP setting	Description
Auto plotting	Enabled	Automatic plotting of a spectra after each experiment. Ensure that this box is on the first time AutoTest is run to produce a plot, but it can be off for subsequent normal maintenance mode.

Table 2 Output configuration (continued)

Option	Typical and ATP setting	Description
Print parameters	Enabled	Automatic plotting of a separate page includes a parameter set, pulse sequence, and descriptive text (only if Auto Plotting is selected). Ensure that this box is "on" the first time AutoTest is run to produce a plot, but it can be "off" for subsequent normal maintenance mode.
Graph history	Enabled	Automatic plotting of history graphs at the end of the AutoTest run (only if Auto Plotting is selected).
Process during acq.	Enabled	Automatic (wnt) processing and display of spectra after each FID.
Plot to file (Postscript printer only)	Disabled	Automatic plots are stored as Postscript files in ~/vnmrsys/autotest/data using the name defined by the global <code>at_currenttest</code> followed by <code>.ps</code> . The system printer should be set to a postscript printer prior to starting AutoTest.

These check boxes reflect the setting (y or n) of global parameters `at_plotauto`, `at_plotparams`, `at_graphs`, and `at_wntproc` respectively. The global parameter `at_plot_to_file` is set to "y" for postscript plots and "n" for hardcopy plots. These global parameters, as well as those mentioned for test selection, are updated when the **Begin Test** button is selected.

NOTE

If you select **Auto Plotting**, **Print Parameters**, or **Graph History**, ensure that enough paper is available for the printer or plotter. If Postscript plotting is selected, ensure that the printer is set to a Postscript printer. (Use the Main Menu selection **File > Printers**.)

Enable tests

The check boxes on the AutoTest window set whether VT, ^{13}C , lock, and gradient tests are enabled or not. These check boxes should be set appropriately before selecting the **Begin Test** button. If these are not selected, any tests involving them are skipped in an AutoTest run, even if the tests are selected in a test package or individually.

The normal use of these check boxes is to allow skipping of some tests when the **All Standard Tests** package is used. The boxes reflect the current state of the options when the AutoTest program starts (these boxes reflect the user global parameter values, y or n, of `at_vttest`, `at_13Ctests`, `at_locktests`, and `at_gradtests`, respectively).

Table 3 Options in the AutoTest Window

Option	Typical and ATP Setting	Description
Enable VT test	Enabled	Tests the VT to check if the hardware is present and the probe supports VT operations.
Enable ^{13}C tests	Enabled	Tests ^{13}C related instrument and probe performance.
Enable lock tests	Enabled	Tests lock functions and stability.
Enable gradient tests	Enabled	Test gradient hardware and probe.
Repeat until stopped	Disabled	Automatic repeating of AutoTest (until manually aborted).

In most cases, a single AutoTest run is appropriate and the **Repeat Until Stopped** check box should not be set. For troubleshooting, or to establish a statistically valid database, however, you might want to run the test overnight. When the **Repeat Until Stopped** check box is selected, the automatic processing and display of data after each FID option is disabled after the first pass (because it is unlikely that the user is present to view the data). In addition, all function tests that do not produce a numerical result are skipped. This check box sets the global parameter `at_cycletest`.

Select test packages

In the center of the AutoTest window there is a column of check boxes that is used to select the test package.

The first check box is **All Standard Tests**. This option sets AutoTest to make a full run, checking out all the relevant hardware enabled by the upper check boxes.

NOTE

All Standard Tests must be the first test performed (with all options enabled) because it determines calibrations for any of the specific tests and is required for the NMR Spectrometer ATP.

Click the check boxes next to each output and the test options that are to be used or click the check box next to **All Standard Tests** or **All Standard Tests** (Power Limited) to select the standard test array.

The list of test spans two pages. Click the arrows to move between pages.

Once the full AutoTest is run, any single test may be run and it can use the calibrations stored in the standard parameter set or in the global variables updated by specific tests. Thus, a single test can be run without doing any calibrations, but its accuracy is dependent on the last calibration performed.

As AutoTest experiments are completed, relevant calibrations are stored in special global parameters (see [Table 4](#) on page 49). In addition, the standard.par parameter set (stored in autotestdir+/parameters) is updated whenever the tof, pw90, T1, linewidth, or gain values are determined. These are determined in the first few experiments of the **All Standard Tests** run, or if specific tests are requested.

Thus, if only a single specific test is requested, the standard.par parameters should have appropriate values so that a full auto calibration of all parameters is unnecessary.

The relevant global parameters include ^{13}C -related parameters, gradient calibrations, etc. The macro /vnmr/autotest/macolib/ATglobal creates these variables (if not already present) and you can read it to get an idea of the global variables used by AutoTest. These global parameters are displayed in the AutoTest Settings popups (see [“AutoTest Settings”](#) on page 26).

Below the **All Standard Tests** check box are several check boxes for different packages of tests. By selecting one or more of these, you can choose the tests performed. Selections may be unselected by clicking the check box once again.

If there is not enough room within the display to show all choices, a pair of arrow symbols appears at the top of the list. Select the right arrow to update the display with the next list of choices. Select the left arrow to return to the previous page of choices.

Run AutoTest

A full AutoTest run including all available options must be run before any single test or partial set of tests is specified. Most of these tests rely on calibrations that are performed as part of the full AutoTest run. This option is specified by the **All Standard Tests** or **All Standard Tests (Power Limited)** check boxes.

Do either of the following:

- If you have previously configured all the required experiments and do not want to make any changes, click the **Begin Test** button to start AutoTest.
- If you have not previously configured all the required experiments or want to change the list of experiments see [“Configure AutoTest”](#) on page 38 .

After selecting the **Begin Test** button, AutoTest begins. The total time for the test(s) will depend on the test(s) specified, on plotting, and the CPU speed.

As AutoTest runs, the FIDs are stored in the data directory, and the results from the tests are stored in the history directory.

Standard Tests Performed by AutoTest

This section lists the standard tests performed by AutoTest.

Automated console tests

- 90° and 30° pulse stability
- 1 μ s amplitude and turn on stability
- Pulse turn on time
- DANTE turn on test
- Quadrature image: 1 scan and 4 scans
- Quadrature phase selection: 0, 90, 180, and 270 degrees
- Phase stability test (13° test)
- Phase switch/settling time
- Attenuator test full and reduced power
- Attenuator test for channel 2 as ^{13}C
- Modulator linearity

- Small-angle phase shifting 0–360 degrees
- High-band amplifier compression
- Low-band amplifier compression
- Temperature homogeneity and rise in decoupler heating test
- Temperature increase in spinlock test
- Temperature jump test
- Sensitivity for 90° and 10° pulses
- AutoGain result for 90° pulse
- Receiver gain
- Folded noise reduction with large spectral width
- Benefit of oversampling at normal gain
- Benefit of oversampling at normal gain + 12dB
- Signal-to-noise as function of gain
- Spectral purity (“glitch test”)
- Lock power test correlation coefficient
- Lock gain test correlation coefficient

Automated tests with shaped RF

- Gaussian 90° stability
- Gaussian phase stability test
- Phase-ramped Gaussian phase stability test
- Shaped pulse accuracy- gaussian excitation profile
- RF amplitude predictability using a gaussian pulse at variable power
- RF amplitude predictability using a gaussian, rectangular and eburp-1 pulses
- Amplitude scaling of shaped pulses using a gaussian pulse
- RF excitation predictability using a variety of shaped pulses

Automated decoupling performance tests

- ¹³C phase modulation decoupling profiles
- ¹³C GARP decoupling profile
- ¹³C WALTZ-16 decoupling profile
- ¹³C XY-32 decoupling profile

- ^{13}C MLEV-16 decoupling profile
- ^{13}C adiabatic decoupling profiles (if waveform generator present on decoupling channel)
- ^{13}C STUD decoupling profile
- ^{13}C WURST decoupling profile
- Sample heating during ^{13}C broadband decoupling

Automated 90 degree pulse width calibrations (PW90)

- ^1H 90° pulse width calibrations at high and reduced power
- ^{13}C 90° pulse width calibrations at high and reduced power

RF homogeneity tests

- ^1H rf homogeneity test
- ^{13}C rf homogeneity test

Gradient calibrations and performance tests

- Gradient level for 10 G/cm along the following:
 - Z axis for all gradient probes
 - X axis for triax probes
 - Y axis for triax probes
 - Cancellation following a gradient
- Gradient echo stability for the following:
 - Z axis at 30 G/cm
 - X axis at 10 G/cm
 - Y axis at 10 G/cm
 - Z axis at 10 G/cm
- Gradient recovery stability for the following:
 - X axis at 10 G/cm
 - Y axis at 10 G/cm
 - Z axis at 10 G/cm
- Gradient recovery (rectangular and shaped gradients) for the following:

X axis at ± 10 G/cm

Y axis at ± 10 G/cm


Z axis at ± 10 G/cm

- Cancellation after gradient
- CPMG T_2 result for the following

Gradient level = 10 G/cm

Without gradients

1% gradient mismatch



4 AutoTest Administration

Saving Data and FID Files from Previous Runs	48
Creating Probe-Specific Files	48
AutoTest Directory Structure	49
AutoTest Macros	54



Saving Data and FID Files from Previous Runs

As AutoTest executes, data and FID files are written into the history and data directories, which are located in the autotest directory. The autotest directory is usually located in the directory `vnmrsys` of a user's home directory. The contents of the data directory are progressively overwritten as AutoTest continues.

Before starting a new AutoTest run, do the following to save the data from a previous run:

- 1 Open a LINUX/UNIX window and enter
`cd ~/vnmrsys/autotest.`
- 2 Change the name of the history directory by entering, for example, `mv history history.old.`
- 3 Change the name of the data directory by entering, for example, `mv data data.old.`

Creating Probe-Specific Files

If you run AutoTest with different probes, you should keep separate autotest directories. Use the following steps to create probe-specific files.

- 1 After you have run AutoTest using a specific probe, change the name of the autotest directory by using the `mv` command, for example:

```
cd ~/vnmrsys
mv autotest autotest_probe_1
```

Where *probe_1* is the name of the probe that was tested, for example, `5mmTriplePFG` or `5mmID`.

Any new AutoTest run automatically creates a new autotest directory in the user's `vnmrsys` directory. The only file that needs to be updated would be `~/vnmrsys/autotest/parameters/standard.par`. This can either be copied from the saved autotest file or the parameter set may be retrieved using `rt` or `rtp`, the parameters updated and then saved, replacing the `standard.par` file. This should be safe for most parameters, but there are several parameters dealing with gradients and indirect detection that must also be checked. It is safest to do an All Standard Tests run the first time a new probe is used.

Once a calibrated `standard.par` parameter set is present, autotest directories may be renamed whenever a probe is changed. In this way, history files may be kept specific to a probe.

- 2 To change the file name back to `probe_1` (or the name you have chosen), enter, for example:

```
cd ~vnmrsys
mv autotest autotest_probe_2
```

Where `probe_2` is the name you have chosen for the probe last tested.

```
mv autotest_probe_1 autotest
```

Where `probe_1` is the name you have chosen for the probe you now want to test.

If you need to repeat any individual test, you can do so by recalling the appropriate FID from the `data` directory. The experiment can then be started with the `go` command without overwriting the previous data. Or the test may be selected from the **Test Library** after using the autotest macro or a menu calling this macro.

AutoTest Directory Structure

AutoTest uses the `~/vnmrsys/autotest` directories listed in [Table 4](#).

Table 4 AutoTest directories

Directory	Contents
<code>data</code>	FIDs from the recent AutoTest run(s)
<code>data.failed</code>	FIDs from any failed Auto Test experiments
<code>history</code>	History files for the various tests
<code>reports</code>	Copies of the report generated each time AutoTest is run
<code>parameters</code>	Parameter files-default entry is <code>standard.par</code>
<code>texts</code>	Copies of the text files attached to the AutoTest experiments
<code>atdb</code>	AutoTest database

data directory

The `~/vnmrsys/autotest/data` directory contains FIDs collected in previous AutoTest experiments. As each experiment finishes, the macro specified by the `wexp` parameter executes, and as part of that macro, a `svf` command is performed that saves the FID under a file name specified by the parameter `at_currenttest` (if it contains a name). The macro first removes any file by the same name (the results of the test from the last time it was run) and then executes `svf`. Thus, the data directory may contain FIDs obtained during different AutoTest runs if those runs were not full runs.

Data files stored in the data directory can be recalled by normal VnmrJ commands such as `rt`. The data may then be transformed and displayed. The `wexp` parameter will contain the name of the macro normally used for data processing, so that the `wexp` command can be used to duplicate the actions normally done in an automatic manner.

The data directory also stores several text files that show calibration results and postscript plot files that are stored when a Postscript printer is used and `at_plot_to_file="Y"`. The plotting action is performed by the macro `ATpage`.

NOTE

If the file `~/vnmrsys/autotest/atdb/at_selected_tests` is empty, only processing and no further acquisition is done. If the file `~/vnmrsys/autotest/atdb/at_cycled_tests` is not empty, those tests may start. Therefore, clear the contents of `at_selected_tests` and `at_cycled_tests` before manually executing the macro.

This result is normally the case if the last AutoTest run came to a normal completion. However, if the last AutoTest run was aborted and no new entry into the AutoTest Program was done, this file will contain entries and an acquisition may start up following the `wexp` command. If so, just abort the acquisition.

data.failed directory

The `~/vnmrsys/autotest/data.failed` directory contains any data from any failed experiment. Failure results when a calculated result falls outside limits defined in the `~/vnmrsys/autotest/atdb/at_spec_table` file. Agilent specifications are indicated in the `~/vnmrsys/autotest/atdb/at_spec_table` file. Users can modify this file by supplying upper and lower limits. Any user-modified `at_spec_table` file should be backed up

outside `~/vnmrsys/autotest`, since this file can be deleted later.

parameters directory

The `~/vnmrsys/autotest/parameters` directory contains any parameter set used by AutoTest macros, including any put there by the user. Normally, only `standard.par` is present. This parameter set has all parameters necessary for the AutoTest macros. Some parameters are only displayed when certain variables are nonzero, or "y" if a string parameter; however, these parameters are printed and displayed if used in an experiment. The AutoTest macro `ATrtp` is used to recall a parameter set from this directory.

reports directory

The `~/vnmrsys/autotest/reports` directory contains text files from previous AutoTest runs, by date. Each run produces a report, whether plotting is requested or not. The report file for a currently proceeding AutoTest run is `~/vnmrsys/autotest/REPORT`. At the end of an AutoTest run, this file is copied to the reports directory under a title that includes the date and time. If AutoTest is repeated, a new report is automatically written out for each complete AutoTest run. The existing `~/vnmrsys/autotest/REPORT` file is renamed as `~/vnmrsys/autotest/LASTREPORT` whenever an AutoTest run begins. Similar actions are executed for the `atrecord_report`. This directory also stores any `FAILREPORTs` with appropriate date-stamps.

texts directory

The `~/vnmrsys/autotest/texts` directory contains mainly text files that are printed on some spectral plots and most parameter set printouts. These files explain the purpose of the test.

history directory

The `~/vnmrsys/autotest/history` directory contains text files that record the values determined in AutoTest runs. They are generated automatically by the `ATrecord` macro, which is used in any AutoTest macro that obtains a numerical result from an NMR experiment. Each result is written on a new line and is date-stamped. Tests that have an Agilent specification listed in the manual *Acceptance Test*

Procedures will be denoted as having passed or failed.

If the history file has more than one result per line, any one failure will cause a fail result for the whole line. When the history file is viewed using the History display (after using the macro `autotest`), failure is indicated by a red data point in graphical output and a colored entry in the text output.

If a user writes a new AutoTest macro including the `ATrecord` macro, the `at_spec_table` must be updated for the history files to be displayed. In addition, the new macro must be listed within the `at_tests_file`.

atdb directory

The `~/vnmrsys/autotest/atdb` directory contains mainly the following text files used by the Auto Test program to create the AutoTest interface:

at_tests_file

The `at_tests_file` file defines all the tests that AutoTest can perform. Tests are specified by a macro name and description. Normally, these are grouped and separated by a line starting with `Label`. The word following it will be displayed as a heading for a group of tests. The test descriptions are displayed in the Test Library display. The macro names are not shown in the display, just check boxes next to the test description.

New tests may be added to the `at_tests_file` by specifying a group title (use the `Label` keyword as the first word on the line, followed by a descriptive phrase). Specify a macro name and then a test description, one per line.

at_groups_file

The `at_groups_file` file defines test packages that have been assembled for convenience. Each package has a line that gives a description (in double quotes), followed by a list of macros to be used in the order of acquisition. There are no restrictions on the placement of these macros in the text file; only the order matters. When the next double-quoted entry appears, a new group is set.

The AutoTest interface display shows these packages as check box entries in the **Configuration** display. Selection of one or more of these causes their execution in the order of selection, once the **Begin Test** button is selected. When this

happens, the `at_selected_tests` file is fixed. Selection of the **All Standard Tests** check box disables any other selections that will be done as part of the **All Standard Tests** run.

Users may add new packages to the **Configuration** display list by adding appropriate lines to the `at_groups_file` in the same format.

Any macro specified within a group must be defined in the `at_tests_file`.

at_selected_tests

The `at_selected_tests` file contains the names of the macros to be run as part of the AutoTest procedure and is fixed at the time the **Begin Test** button is selected. The format is one line per macro with each line containing the name of a macro, in the order of acquisition.

As AutoTest proceeds, each line is deleted as the specified macro finishes its activity. Thus, completion of the AutoTest run is defined as when this file is empty. The single exception is the case of automatic repeating of AutoTest, as specified by the **Repeat Until Stopped** check box in the **Configuration** display and as indicated by the value of the global parameter `at_cycletest ("y")`. In this case, at the completion of the AutoTest run, the contents of the file `at_cycled_tests` are copied into `at_selected_tests` and the process then continues.

at_cycled_tests

The `at_cycled_tests` file is updated when the **Begin Test** button is selected. If the **Repeat Until Stopped** check box is selected in the **Configuration** display, the global parameter `at_cycletests` is set to "y" and the file `at_selected_tests` is copied to `at_cycled_tests`. If no test cycling is requested, this file is emptied.

at_spec_table

The `at_spec_table` file is written out when the `~/vnmrsys/autotest` directory is created and is spectrometer-dependent. The appropriate file is copied from the directory `/vnmr/autotest`, depending on the Spectrometer frequency. It contains a list of macros used in AutoTest for producing entries in the history files. For each macro the following is specified:

- The history file affected by the macro.

- The column number (not counting date) containing the result.
- The lower limit for the result.
- The upper limit for the result.
- A text description of the result. This text description is used for the graphical displays and plots. A comment line above each macro serves to describe the test.

All results specified in the manual *Acceptance Test Procedures* have upper and lower limits specified numerically in this file. Those not having Agilent specifications have asterisks (*) as entries for upper and lower limits and these results will have no indication of pass or fail in their history files, or colored indication of failure in the graphical displays of the history files.

Users may wish to set their own upper and lower limits for many, if not all, of the results. They may do so by replacing the asterisks with numbers. Of course, this should only be done after a good statistical base is obtained, such as more than 20 complete AutoTest runs. Once this base is obtained, the numbers put into the `at_spec_table` file should have a reasonable margin of error built in.

It is a good idea to make a copy of the `at_specs_table` file prior to changing it, as well as the modified file; because, deletion or renaming of the `autotest` directory will result in a default `at_spec_table` being copied from `/vnmr/autotest/atdb`.

AutoTest Macros

To help users who want to add or modify tests, this section describes some of the macros used by AutoTest. These macros are in `/vnmr/autotest/maclib`.

ATglobal macro

The `ATglobal` macro is run when the AutoTest program begins. The macro checks for the existence of `autotest` parameters in the user file `~/vnmrsys/global`. These parameters are used to store calibrations and results that are used by AutoTest macros. If the parameters are not present, `ATglobal` creates them. Otherwise, the parameters are left unchanged. In VNMRJ, the **Tools > Standard Calibration Experiments > AutoTest Settings** drop-down menu option permits

convenient viewing (and entry of) AutoTest global parameters. A partial list of these parameters is given in Table 5.

Table 5 Selected Parameters Created by ATglobal

Parameter	Contains
at_currenttest	Name under which the FID is stored
Autotestdir	Full path of the autotest directory
at_user	Name of the user running autotest (printed in the report)
at_coilsize	length (in mm) of active window in coil (typically 16 or 18 mm)
at_consoletype	Name of the console entered in AutoTest window
at_consolesn	Number of the console entered in AutoTest window
at_probetype	Name entered for probe used in AutoTest window
at_wntproc	y or n (for processing/display after each FID)
at_cycletest	y or n (for automatic repeating of AutoTest)
at_printparams	y or n (for parameter list/pulse sequence printouts)
at_plotauto	y or n (for automatic plotting)
at_plot_to_file	y for (for Postscript plotting to data directory)
at_graphhist	y or n (for history graphs plotting)
at_locktests	y or n (for lock power/gain tests)
at_max_pwx1vl	Maximum permitted ^{13}C power level
at_max_pwx2lvl	Maximum permitted ^{15}N power level
at_max_dpwr	Maximum permitted ^{13}C decoupling power
at_max_dpwr2	Maximum permitted ^{15}N decoupling power
at_T1	Value of last determined T_1
at_gain	Value of gain determined by auto gain
at_tof	Value of tof for water
at_fsq	Value of fsq parameter
at_dsp	Current value of dsp at the start of run
at_ampl_compr	Value of the high-band amplifier compression

Table 5 Selected Parameters Created by ATglobal (continued)

Parameter	Contains
at_LBampl_compr	Value of the low-band amplifier compression
at_LBampl_compr 10us_c	¹³ C amp compression at at_pwx90lvl_10 sec_c
at_decHeating	Temperature increase from decoupling
at_linewidth	Linewidth of water resonance
at_pw90	90° pw at power specified in the AutoTest display
at_tpwr	Power specified in the AutoTest display
at_pwx90c	¹³ C pw90 at power at_pwx90lvlc
at_pwx90lvlc	¹³ C power level for ~15 µs pwx90 at power <= to at_max_pwxlvl
at_pw90Lowpower	90° pw at reduced power
at_pwx90Lowpowerc	¹³ C pw90 at power at_pwx90Lowpowerlvlc
at_tpwrLowpower	Power level at reduced power
at_pw90_ch2	90° pw on channel 2
at_pwx90	¹³ C pw90 determined at at_pwx90lvl
at_pwx90lvl	¹³ C power level for approximately 15 µs pwx90
at_pwx90Lowpower	¹³ C pw90 at reduced power
at_pwx90Lowpowerlvl	¹³ C power level at reduced power
at_pwx90_10us_c	¹³ C pw90 at power level at_pwx90lvl_10 sec_c
at_pwx90lvl_10us_c	¹³ C power level for ~10 µs pwx90 at power<= to at_max_pwxlvl (typically at 800 MHz)
at_pwx90Lowpower_10us_c	¹³ C pw90 at at_pwx90Lowpowerlvl_10 sec_c
at_pwx90Lowpowerlvl_10us_c	¹³ C power level for at_pwx90Lowpower_10 sec_c
at_vttest	y or n (for VT test)
at_temp	Current temperature
at_vttype	Current value of the global parameter vttype
at_tempcontrol	Value reflects usage of temp tcl/tk panel
at_gradtests	y or n (for gradient tests)

Table 5 Selected Parameters Created by ATglobal (continued)

Parameter	Contains
at_pfgon	Current value of pfgon
at_gmap	y or n (for gradient tests)
at_gzcal	Value of G/cm per DAC unit for z-axis gradient
at_gxcal	Value of G/cm per DAC unit for x-axis gradient
at_gyca1	Value of G/cm per DAC unit for y-axis gradient

ATstart macro

The ATstart macro is run after the **Begin Test** button is selected in the **Configuration** display. The macro sets the global parameters to reflect the current state of the hardware and aborts under certain circumstances, such as if requested tests are not compatible with the current hardware settings.

Messages are displayed indicating the source of the problem. The reports are initialized with relevant information and the ATnext macro is executed.

ATnext macro

The ATnext macro checks the at_selected_tests file and copies the first entry into the global parameter at_cur_smacro, deletes the top line in at_selected_tests and executes the macro specified by at_cur_smacro. If the at_selected_tests file is empty, ATnext either finishes the AutoTest run or calls the ATrestart macro which copies the at_cycled_tests file to at_selected_tests, permitting repeated AutoTest runs until manually aborted by the user.

ATnext is usually found at the bottom of each macro defining a particular test. This permits the linking of one test to another, in a general fashion.

ATxxx macros

Specific tests usually have the designation of AT, followed by a number or group of letters. Each macro is self-contained, having the ability of setting up parameters, performing acquisition, processing the acquired data, possibly setting up new experiments and processing the data acquired from those experiments, creating plots, parameter printouts,

archiving the raw data, performing statistical analyses of the data, and writing results to history files and reports.

To better illustrate the structure of these macros, [Table 6](#) on page 58 gives the source code for macro AT16, the turn-on test (channel 1). A column of descriptive comments has been added to clarify the statements.

ATcxxx macros

ATcxxx macros use the ^{13}C and ^{15}N power limits.

ATrecord macro

The ATrecord macro is run whenever a numerical result is stored in a history file. It is a general macro that will create the specified history file if it is not present. The macro has a minimum of four arguments: Name of history file, comment line, column header line (name of the result), and value of result. For example, see the above macro near the bottom. The variables \$turnon and \$corrcoef are calculated prior to using the ATrecord macro, and these appear in two columns headed by time and corr_coef. Note the use of trunc. This is necessary to limit the number of decimal places produced.

Up to seven results may be stored in one history file. If this macro is used, be sure to modify the at_spec_file in ~/vnmrsys/autotest/atdb to add the necessary number of lines describing the results and any upper and lower limits desired. If a new history file fails to appear, it is usually because of the failure to update the at_spec_file. A backup copy should be made of the atdb in case of accidental overwriting.

The AutoTest macros can be run as independent macros if a specific test is desired. This is done by either entering the macro in the VnmrJ 3 command line or by selecting the single test from the AutoTest test library panel by using a check box and the **Begin Test** button. Again, if the file at_selected_tests is not empty, the ATnext macro will start a new acquisition.

Table 6 Source code for AT16 macro example

if (\$#=0) then	First time AT16 is run it has no arguments.
ATrtp('standard')	Recalls standard parameter set.

Table 6 Source code for AT16 macro example

text('Pulse Turnon Test')	
at_currenttest='turnon_ch1'	Puts name of test in global variable.
tpwr=tpwr-6 ph	
getlimit('tpwr'):\$max,\$min,\$step	
if (\$step=1) then	
setlimit('pw',100,0,.025)	
array('pw',37,0.1,.025)	Sets up pulse width array.
else	
setlimit('pw',100,0,.0125)	
array('pw',37,0.025,.0125)	
endif	
ss=2	
wnt='ATwft select(celem) aph0 vsadj dssh dtext'	Specifies what to do every FID
wexp='AT16(`PART1`)'	Specifies what to do at end of experiment.
ATcycle	Disables wnt processing if in repeat mode.
au	Begins acquisition and specifies wnt/wexp processing to occur.
write('line3','Pulse Turnon Test ')	
dps	
elseif (\$1='PART1') then	This part executes at end of experiment.
wft	
if (at_plotauto='y') then	
if (at_printparams='y') then	
pap ATplttext	If parameter printout is requested.
pps(120,0,wcmax-120,90)	
ATpage	
endif	
endif	
select(arraydim) aph0	
f peak:\$ht,cr rl(0) sp=-1p wp=2p vsadj dssh	
dtext	
ATreg6	Fits to straight line and displays/plots data.

Table 6 Source code for AT16 macro example

ATpl3:\$turnon,\$corrcoef	Determines turn-on time and correlation coefficients
\$turnon=trunc(\$turnon)	Limits number of decimal places.
\$corrcoef=trunc(1000*\$corrcoef)/1000	
ATrecord('TURNONch1','Pulse Turnon Time (nsec) ',	Writes out results to history file.
'time ', \$turnon, ' corr_coef.', \$corrcoef)	
write('file', autotestdir+'/REPORT',	Writes results to report.
'Pulse Turnon Time : %2.0f	
nsec.-Corr. Coef. = %1.3f ', \$turnon, \$corrcoef)	
if (at_plotauto='y') then	
ATplttext(100,wc2max-5)	
full wc=50 pexpl ATpage	Plots regression fit
endif	
ATsvf	Removes old data set and stores FID under name in at_currenttest.
ATnext	Starts next macro in at_selected_tests file, if present.
endif	Closes elseif part



5 AutoTest - Test Reference

RF Performance Test Descriptions	62
Shaped Pulse Test Descriptions	69
¹³ C Test Descriptions	72
Gradient Tests Descriptions	74
CPMG T ₂	75
¹³ C Power-Limited Pulse Tests	77
¹³ C/ ¹⁵ N power-limited decoupling tests	81
¹⁵ N Power-Limited Pulse and Decoupling Tests	88
2H Pulse and Decoupling Tests	89
Installation Tests for Cryogenic Probes	92
Other Test Descriptions	93
Tests Using Salty Sample	95

All units of the rf system (transmitters, linear modulators, rf attenuators, amplifiers, receivers, and probes) must be in the standard configuration when AutoTest is run. If the system configuration has been changed, it must be returned to the standard configuration before running AutoTest for acceptance testing.

All data is stored, and both plots and statistical analyses are provided as part of the acceptance testing. Plots and statistical analyses are made concurrently with acquisition.



RF Performance Test Descriptions

This section describes the RF performance tests in AutoTest.

Pulse Stability and Sensitivity

Experiment - A single-scan pulse experiment is repeated 20 times and the spectra plotted in a horizontal stack. The average peak amplitude and the rms deviation are measured and reported. This test is run for the following:

- 90° flip pulses
- 30° flip pulses
- 10° flip pulses
- 1 μ s pulses

Purpose of 90° Pulse Stability - Modern experiments require very high pulse reproducibility to minimize cancellation residuals and T_1 noise in 2D experiments. This test checks sensitivity and amplitude reproducibility by comparing a series of spectra obtained with the signal following a single 90° pulse. The statistical analysis produces an rms deviation, in percent of the average peak height.

Purpose of 10° Flip Pulses - The 10° flip data are acquired in the same manner with an additional 12 dB of gain. This is done to ensure enough gain so that the S/N is not dominated by ADC - round off noise (significant at low gain) and that S/N can be compared for different probes (cold vs. warm) or different fields.

Purpose of 30° Pulse Stability - The sinusoidal nature of the excitation profile makes the signal generated following a 90° pulse less sensitive to error than the signal following a much smaller flip angle pulse (the top of a sine wave is broad and changes in amplitude less for small changes in flip angle than for a smaller pulse). A 0° flip angle would have the highest sensitivity to flip angle, but would give no signal, of course. A compromise between the extremes of large signal following a 90° pulse, and no signal following a 0° pulse is to use a 30° pulse. The rms deviation is measured from an array of spectra obtained using 30° pulses.

Purpose of 1 μ s Pulse Stability - This test emphasizes the turn-on characteristics of the pulse. Any instability of the pulse rise should give a corresponding reduction of measured stability. Since the flip-angle is much less than a 90° or 30° pulse, the measured stability may be lower. The

rms deviation is measured from an array of spectra obtained using 1 μ s pulses.

Cancellation test

Experiment - Four single-scan, four 2-scan, and four 4-scan 90° pulse spectra are acquired in which the transmitter phase is held fixed and the receiver is phase-cycled 0, 2, 1, 3. Data are plotted in a horizontal stack with the single-scan spectra on scale. The vertical scale is increased by 100 times and plotted in the same manner. Average residual signal for 2-scan and 4-scan cancellation are determined.

Purpose - Modern experiments (HMQC, HSQC, NOE-difference, and so on) often rely on phase-cycling to achieve desired results. This test compares single-transient response versus two- and four-transient response in which the phase-cycling is set to achieve cancellation. Cancellation performance is often dominated by environmental factors (vibration, for example) and the performance of the lock circuit. Amplitude stability tests (see above) measure signal reproducibility in single scans, while cancellation involves multiple scans, permitting environmental effects to dominate the result.

Phase stability (13 deg. phase error) test

Experiment - The 90° pulse stability test is repeated but uses a 90° pulse-1 ms- 90° pulse train with the carrier positioned 37 Hz off-resonance from the water.

Purpose - Phase stability is essential for high-performance modern experiments. Poor phase stability would produce poorer water suppression and increase T_1 noise in 2D NMR. The most robust tests of phase stability are solids tune-up sequences used for verifying performance for line-narrowing sequences, such as WAHUHA or MREV-8, because these sequences are fairly independent of amplitude stability.

This test is the “13° test” in which two 90° pulses separated by

1 ms are applied with the transmitter placed 37 Hz off-resonance. The resulting NMR response stability is a product of both rf amplitude stability and phase stability because variations in phase between the pulses induce an amplitude change. The observed amplitude error should be divided by a factor of 7.1 to obtain a measure of phase error, in degrees.

Pulse turn-on time

Experiment - Single-scan experiments are taken in which the pulse is varied in minimum pulse-width steps at low enough power so that the response is linear. The response is fitted to a straight line and the turn-on time is determined.

Purpose - The quality of modern rf is good enough that examination of pulse shapes using an Oscilloscope is not as informative as well-designed and executed NMR tests. The turn-on and turn-off characteristics of a very short pulse are properties that can be measured sensitively by NMR.

The turn-on test measures the amplitude of a signal following a short variable-length pulse. In the limit of a small flip angle, this dependence is linear. The data are analyzed and least-squares fitted to a straight line. The intercept is the pulse turn-on time.

Attenuator linearity

Experiment - For a small flip angle pulse, the rf coarse power is varied from maximum to minimum in a single-scan mode. The data are plotted in a horizontal stack to facilitate visual inspection. The data are fitted to a linear regression and plotted in phased mode to show any phase change as a function of power.

Purpose - Overall power control is accomplished using PIN diode-controlled rf attenuators. These attenuators are precision devices that should have negligible phase change throughout their full range. The amplitude response should also be logarithmic. A log regression analysis should show the extent of fit to the ideal. The phase change as a function of power is examined.

This test does not permit a full assessment of the cause of the phase error, because the amplifier might be in compression at the maximum power output.

Attenuator linearity at reduced power

Experiment - The attenuator linearity test is repeated but with output of the transmitter reduced by the linear modulator. This is done to isolate the effect of the rf amplifier.

Purpose - The attenuator linearity test is performed, but with reduced power input to the attenuator (using the linear

modulator to reduce the output power from the transmitter).

Linear modulator linearity tests

Experiment - With the coarse rf amplitude set at a value 23 dB down from maximum, the rf power is varied using the linear modulator. The linear modulator is used for fine power control and shaped rf excitation. The rf amplitude is varied, over a range of 60 dB, in 100 equally-spaced steps over the whole range. This should produce a linear ramp of signal response following a small flip-angle pulse. Spectra are plotted in a horizontal stack of spectra in phased mode with the highest signal full scale. The width of the plotted region around the water is set narrow enough to clearly show the base of the water peak. The data are fitted to a straight line using a linear regression analysis and plotted.

Purpose - Further power control is possible using the linear modulator present on the NMR transmitter board. This test produces a series of experiments in which the pulse power is changed over the full range of the modulator. The linear nature is tested by a linear least-squares fit of the data.

Predictable power control is essential for delivering accurate shaped pulses and for precise power level control in Hartmann-Hahn experiments in both solids and liquids.

Linear modulator linearity tests with attenuators set to low power

Experiment - The linear modulator linearity test is repeated, with the coarse rf amplitude control set for very low power. The pulse width is increased correspondingly to obtain comparable signal-to-noise as in the linear modulator linearity test. The data are fitted to a straight line using a linear regression analysis and plotted.

Pulse shape test-DANTE

Experiment - The rf amplitude is set for a 20 μ s 90° pulse, and the result is compared to that for a single-scan spectrum:

- 10 pulses, 2 μ s each
- 20 pulses, 1 μ s each
- 25 pulses, 800 μ s each
- 50 pulses, 400 μ s each
- 100 pulses, 200 μ s each

For all except the first 20 μs pulse, a 1 μs delay is inserted between each pulse. The data are plotted in a horizontal stack to permit comparison of amplitudes. The amplitudes are measured and printed.

Purpose - A DANTE-type test is performed in which the signal response following a 20- μs pulse is measured. This is compared with a series of experiments in which the pulse is increasingly divided into series of pulses spaced by 1 μs . The sum of the pulses is held constant at 20 μs .

If the pulse shape is *ideal* and the total time of the pulse train is short compared to T_2 , the rotation of magnetization should be identical. As the pulse length shortens, any non-ideality of pulse shape is revealed as a reduction in intensity.

Phase switch settling time

Experiment - Parameters p1 and pw are set to the same value

(1 μs) and 30 spectra are acquired using a 2-pulse sequence, with the first and second pulses separated by a delay of 20 μs . The phase of the second pulse is shifted 180° from the first pulse at a variable time prior to the second pulse. A single-pulse spectrum is also acquired. The arrayed 2-pulse spectra and the single-pulse spectrum are plotted, with the single-pulse spectrum last. This last spectrum serves as a reference. The phase shift should be accomplished in 100 ns or less.

Purpose - This test exercises the phase-shift hardware by finding the time needed to perform a 180° phase shift. The pulse sequence is a version of *jump-and-return* in which two

1- μs pulses are executed just 20 μs apart. Ideally, because the second pulse has a 180° phase shift with respect to the first pulse, there should be no excitation. By varying the time before the second pulse, at which the phase shift is done, from 0 to 20 μs , an estimate of the phase switch and settling time can be made. The last spectrum is that from just a single, 1- μs pulse, and serves as a reference. This phase shift should be accomplished in 100 ns or less.

RF homogeneity

^1H RF Homogeneity Experiment - One hundred experiments are run in which the pulse width is incremented. The

spectra are plotted in a horizontal stack in phased mode and sufficiently expanded so that the base of the water can be examined using the same phase settings for each spectrum.

13C RF Homogeneity Experiment - In the pulse sequence `delay-pw90(1H)-delay(1/2JCH)-pw(13C)` the `pw(13C)` is varied from 0 to a flip angle greater than 900° while observing the ^{13}C -coupled protons. The 0-flip-angle spectrum is adjusted to full scale and the data expanded to show only the ^{13}C -bound protons side-by-side to permit measurement of X-coil rf homogeneity. The results are plotted and displayed in magnitude mode showing one of the lines of the methanol doublet.

Purpose - This test checks the homogeneity of the rf field strength throughout the active region of the sample. In an ideal case, for nuclei having reasonable T_2 values, the signal generated following a $360^\circ + 0^\circ$ pulse should match that following a 0° pulse. The signal strength as a function of flip angle should be sinusoidal. The amount of drop off is related to the inhomogeneity of the rf field (RF homogeneity results are shimming-dependent. If results for a given probe are outside of specifications, re-shim the sample or use the sample specified in the probe test manual).

High rf homogeneity is important because many important pulse sequences use a large number of pulses. The signal losses accumulate with each pulse such that, in worst cases, the desired signal is lost. Most heteronuclear, indirect detection experiments on large molecules use HSQC pulse sequence components. These contain 6 to 10 ^1H pulses, including 4 to 8 \times nucleus 180° pulses. High rf homogeneity is especially important in these cases.

Receiver test

Experiment - Single scan spectra are collected that span the range of receiver gain and divide that range into at least 25 evenly spaced values of gain, including the highest and lowest gain values. The data are plotted with the highest signal on scale so that the heights can be easily compared.

The results are fitted to a straight line using linear regression, and the fitted data are plotted. Next, the data are normalized and plotted with the water signal held to a constant height so that the noise levels are easily compared (a few mm of noise in the baseline are provided).

The signal-to-noise ratios for the water line in all spectra

are measured.

Purpose - Receiver gain is selectable in a logarithmic manner (in dB). In an ideal case, variation of receiver gain should produce a logarithmic dependence of signal strength. As the gain is lowered, the noise becomes dominated by the noise generated in the ADCs, not in the preamplifier and probe. Regardless of the signal strength, operation in this range of gain will produce poorer signal-to-noise.

Image rejection test

Experiment - Plot the data from the following tests first in a horizontal stack, with the single-scan data on-scale, and then with the vertical scale increased 100 times. Quantitate the average image and center glitch.

- Four single-scan and four 4-scan 90° pulse spectra are acquired in which the carrier frequency is shifted 1000 Hz from the water. The carrier position is not changed during the pulse sequence and acquisition.

Purpose - This test checks for quadrature image. Four single-transient and 4 four-transient responses are collected and compared. A Direct Drive system does not use quadrature detection and will have no images.

Quadrature phase selection

Experiment - Plot the data from the following tests first in a horizontal stack, with the single-scan data on-scale. Quantitate the average intensities.

- With constant receiver phase, acquire multiple spectra with the observed pulse set to a phase of 0, 90, 180, and 270 degrees. Phase the 0 degree spectra to positive absorption. The spectra should show absorption, dispersion (+), absorption (-), and dispersion (-). Calculate the average and standard deviation of the spectra grouped by the phase.

Purpose - This test checks the quadrature phase shift which is used within any phase-cycling experiments.

Shaped Pulse Test Descriptions

This section describes the shaped pulse tests for AutoTest.

Gaussian-shaped pulse excitation

Experiment - A gaussian-shaped pulse, with excitation bandwidth at 50% amplitude about 200 Hz, is applied (for example, a 12 ms, 90° pulse length). Single-scan spectra are taken with the transmitter stepped over the range ± 250 Hz from resonance, in 5-Hz steps.

The data are plotted in a horizontal stack, with the on-resonance spectrum at full scale to illustrate the gaussian shape of excitation. The vertical scale is increased by $\times 10$ and plotted again to show the wings.

Purpose - The most demanding test of shaped pulse accuracy is the ideality of the NMR data following a shaped pulse. This test determines the accuracy of a gaussian pulse by examination of the off-resonance excitation. This is done by repeating the same single-pulse excitation while varying the transmitter position through a wide range.

A stacked array of data should show the magnitude of excitation as a function of offset from the resonance. In the ideal case, this excitation envelope would also be gaussian. Any non-gaussian nature of the pulse, as *delivered to the probe*, would be represented by a convolution of excitation envelopes. For example, if the power were not delivered in a linear manner, producing some rectangular nature, the excitation envelope would have some $\sin x/x$ nature, producing characteristic sinc wiggles. The lack of such non-gaussian behavior is a direct measure of the accuracy with which the hardware can deliver an ideal shape to the nuclei.

Gaussian 90 degree pulse stability

Experiment - The rf 90° pulse stability test is repeated using a gaussian pulse at the same peak power.

The data are plotted in a horizontal stack, with the on-resonance spectrum at full scale to illustrate the gaussian shape of excitation.

Purpose - Modern experiments require very high pulse

reproducibility to minimize cancellation residuals and T_1 noise in 2D experiments. This tests amplitude reproducibility by comparing a series of spectra obtained with the signal following a single 90° pulse. The statistical analysis produces an rms deviation, in percent of the average peak height.

Gaussian 13 degree phase error

Experiment - The rf 13° test is repeated using a gaussian pulse at the same peak power as in the phase stability test.

The data are plotted in a horizontal stack, with the on-resonance spectrum at full scale.

Purpose - The rf 13° test can be done using shaped pulses. In this case, a gaussian pulse is used at high peak power.

Gaussian SLP 13 degree phase error (phase-ramped gaussian pulses)

Experiment - This 13° test is repeated using a phase-ramped gaussian pulses. The rf carrier should be 37 Hz off-resonance from water, but the center of excitation of the gaussian phase-ramped pulses should be 1000 Hz from the carrier. The amplitude of the gaussian pulses is set high enough to exert a 90° pulse on the water.

Purpose - The 13° test can be done using phase-modulated pulses. These type of pulses provide single- or multiple-frequency selective excitation through the use of both amplitude and phase modulation.

Shaped pulse settability

Experiment - Single-pulse, single-scan spectra are collected. The rf power is dropped in eight successive spectra by 3 dB each time and the pulse width increased so that a 90° flip angle is maintained. The spectra are plotted in a horizontal stack for easy amplitude comparison.

Purpose - An rf attenuator should permit accurate power control. In this case, the pulse length is repeatedly incremented while appropriately reducing power levels. The NMR response should be identical.

Shaped pulse test - rectangular, gaussian, and EBURP-1

Experiment - Single-scan, one-pulse excitation spectra are collected using rectangular, gaussian, and EBURP-1 pulses at

the same peak amplitude (note the power value and pulse lengths). Constant peak amplitude is maintained; therefore, pulse width ratios of 1.0:2.4:16.0 for the rectangular: gaussian: EBURP-1 pulses, respectively, are used to obtain the same flip angle. Spectra are plotted side-by-side in absolute intensity mode at full vertical scale.

At any constant power, the 90° pulse lengths should reflect their theoretical ratios. Here, the pulse lengths are set in a ratio of 1:2.4:16. The resulting NMR responses should be identical in amplitude.

Shaped pulse test-constant bandwidth for a variety of shapes

Experiment - Single-scan, one-pulse excitation spectra are collected using a variety of shapes that are automatically calculated using Pbox, based on a single pulse calibration using a rectangular pulse, for a constant 4000 Hz bandwidth. The shaped pulses have different peak amplitudes and pulse widths (note the power value and pulse lengths). Spectra are plotted side-by-side in absolute intensity mode at full vertical scale. The resulting NMR responses should be identical in amplitude.

Shaped pulse scalability

Experiment - A small flip-angle gaussian pulse is used for a single-transient, one-pulse spectrum. The linear modulator is used to scale down the amplitude of the pulse in 100 steps over a range of 60 dB. Plot widths are set small enough to show the base of the water and plot all spectra in a horizontal stack in phased mode with the maximum signal spectrum at full scale.

Purpose - The linear nature of the system is graphically tested by measuring NMR response when the amplitude is under full control, both by the rf attenuator and by the linear modulator.

13C Test Descriptions

X-coil rf homogeneity can be determined using an indirect detection pulse sequence. Sensitivity in many indirect detection experiments is markedly affected by X-coil performance because of the large number of 180° pulses used.

X-decoupling is tested for various modulation schemes at constant amplitude (WALTZ-16, GARP-1, and so forth.) as well as more powerful adiabatic pulse techniques. Efficiency is measured by varying the 13C decoupling frequency while observing the proton spectrum under broadband decoupling.

13C 90 degree pulse width calibration

The power level for a 90° flip of approximately 15 μs on the X-coil of the probe is determined. Amplifier compression is determined by lowering power by 12 dB and redetermining the 90° pulse width. Both results are reported.

X-Phase modulation decoupling profiles

13C power level is reduced 20 dB from the level used to obtain a 15 μs 90° (approximately 1.8 kHz), and the 13C decoupling efficiency is determined for the following phase-modulated, constant-amplitude broadband decoupling sequences:

- WALTZ-16
- GARP-1
- XY-32
- MLEV-16

The 13C decoupling frequency is varied over a range of ± 80 ppm in a series of single-scan, proton-observe experiments. Only the 13C-bound protons are shown in the expanded spectrum, which is plotted with the spectra side-by-side in absolute intensity mode to illustrate decoupling efficiency.

13C-Adiabatic decoupling profiles

The decoupling profile experiment is repeated with the following adiabatic decoupling schemes:

- STUD modulation
- WURST modulation

Decoupler Heating

The same test as in the variable temperature test is performed but this time using a 75 ms ^{13}C decoupling period prior to acquisition within a total recycle time of 1.5 seconds, including acquisition. One-hundred, single-scan spectra are collected with ^{13}C decoupling followed by 100 identical spectra with no decoupling. The spectra are plotted in a stacked manner to permit examination of the rate of change of temperature, the homogeneity of temperature, and the length of time necessary to reach equilibrium. The rf field strength must be sufficient to decouple over a 160 ppm range using garp-1.

Gradient Tests Descriptions

This section describes the gradient tests in AutoTest.

Gradient profile

Experiment - A spectrum with a 100 kHz spectral width is acquired using a gradient echo (collect echo during a Z-axis gradient). This acquisition is repeated for both positive and negative gradients that are sufficient to spread the pattern greater than 50 kHz at the base. Gradient strength and duration as well as the size of the active length of the coil are noted. The experiment is repeated for both the X-axis and Y-axis gradients, if available.

Purpose - This test uses pulsed field gradients (PFGs) to quantitate the gradient field strength. The width of the pattern is directly proportional to the gradient strength. The width at 20% of maximum is used to calculate the gradient strength. Both positive and negative gradients are used. This should be done for all orthogonal axes available.

Field recovery stability

Experiment - The 90° pulse stability test is performed but, preceding the rf pulse, is a 1 ms, 30 G/cm Z-axis gradient pulse, which is then followed by a 100 μs field stabilization delay. This test is repeated with a 10 G/cm gradient pulse. If X and Y gradients are available, the test is repeated with a 10 G/cm gradient pulse for each gradient.

Purpose - A gradient is applied prior to a measuring pulse. The stability of the signal response is used to measure the ability and reproducibility of the system to recover from the gradient pulse.

Field recovery

Experiment - The gradient 90° pulse stability test is repeated with the field recovery delay varied from 0 to 1000 μs, using a positive 10 G/cm Z-axis gradient pulse. The spectra are plotted in a horizontal stack with the 1000 s data at full scale. The test is repeated using a negative 10 G/cm gradient and, if available, for both the X and Y gradients. Rectangular and half sine gradients are used. Recovery is defined as the time it takes to recover to 95% or more of the amplitude at 1000 μs.

Purpose - A gradient is applied prior to a measuring pulse and the time before the pulse is varied. The rate of recovery determines how soon a pulse may be applied.

Gradient echo stability

Experiment - The 90° pulse stability test is run, this time with a positive gradient for 1 ms, a 500 μ s delay, and a negative gradient for 1 ms following the rf pulse. The following gradient strengths and axes are used.

- 30 G/cm-Z axis only
- 10 G/cm-Z axis and, if present, Y and Z axes

Purpose - Following a single pulse, a pair of oppositely-signed gradients is applied. The stability of the resulting refocused signal measures any instability in the gradient amplitude, as well as the accuracy of the gradient level control. This should be done for all the orthogonal axes available.

Gradient effect on cancellation test

Experiment - Four 1-scan, four 2-scan, and four 4-scan 90° pulse spectra are acquired, with a 10 G/cm Z-axis gradient pulse 100 μ s prior to the rf pulse and the transmitter phase held constant while the receiver is phase-cycled 0, 2, 1, 3. The spectra are plotted in a horizontal stack with the single-scan spectra on scale. The vertical scale is increased by 100 times and plotted in the same manner. The average residual signal for the 4-scan cancellation is quantitated and the results plotted and analyzed as above.

Purpose - The cancellation test is done with a gradient pulse applied 100 μ s before the rf pulse. If the lock circuitry and field recovery characteristics are favorable, no deterioration in cancellation efficiency should be noted.

CPMG T_2

Experiment - A CPMG T_2 experiment is performed with 2 ms between 180° pulses for total echo pulse trains ranging from a few milliseconds out to at least $2 \cdot T_1$. This experiment is repeated for the case where 500 μ s 10 G/cm rectangular pulses are placed around the 180° pulses in each echo, as well as the case in which no gradients are used. The values for T_2 are reported for all cases. The experiment is repeated

for the case of a 1% mismatch in gradient amplitude.

Purpose - Following a single pulse, pairs of same-signed gradients are applied within the echoes of a CPMG T_2 pulse echo train. The measured T_2 gives a measure of any instability in the gradient amplitudes as well as the accuracy of the gradient level control. This measurement is compared to an identical experiment in which the gradient amplitudes are set to zero or mismatched by 1%. More rapid diffusion at higher temperature will cause the gradient/no gradient comparison to worsen. Therefore, performance over time should be compared for the same temperature. (Measurements done at $\sim 4^\circ$ C have shown no difference in measured T_2 for the gradient/no gradient cases, indicating that diffusion is responsible for the difference in these cases at higher temperatures).

13C Power-Limited Pulse Tests

This section describes the power limited 13C tests in AutoTest.

13C pw90 (15 microsecond) and lowband amplifier compression

Experiment - The user-limited attenuator setting for a 90° flip of approximately 15 μ s on the X-coil of the probe is determined. This setting is limited to an upper limit of the user-supplied value. At this power level the pulse width is varied to obtain a pw90 value. Amplifier compression is determined by lowering the power by 12 dB and redetermining the 90° pulse width. The results are reported and stored in a separate history file from the standard history file (13CPW90c instead of 13CPW90).

13C pw90 (10 microsecond) and lowband amplifier compression

Experiment - The user-limited attenuator setting for a 90° flip of 10 μ s on the X-coil of the probe is determined. This setting is limited to an upper limit of the user-supplied value. Amplifier compression is determined by lowering the power by 12 dB and redetermining the 90° pulse width. The results are reported and stored in a separate history file from the standard history file (13CPW90c_10 sec instead of 13CPW90_10 μ s). This test is only used on systems having a low band amplifier and probe capable of this rf field strength (typically 800 MHz systems).

13C Power for pw90=15.0 microsecond and lowband amplifier compression

Experiment - The optimum user-limited attenuator level for a 90° flip of 15.0 μ s on the X-coil of the probe is determined. The power level is limited to an upper limit of the user-supplied value. At this attenuator level, the fine power modulator is varied to obtain a 90° pulse. Amplifier compression is determined by lowering the power by 12 dB and redetermining the 90° pulse width. The results are reported and stored in a separate history file from the standard history file (13CPW90fc_15 μ s instead of 13CPW90).

¹³C Power for pw90=10.0 microsecond and lowband amplifier compression

Experiment - The optimum user-limited attenuator level for a 90° flip of 10.0 μs on the X-coil of the probe is determined. The power level is limited to an upper limit of the user-supplied value. At this attenuator level the fine power modulator is varied to obtain a 90° pulse. Amplifier compression is determined by lowering power by 12 dB and redetermining the 90° pulse width. The results are reported and stored in a separate history file from the standard history file (13CPW90fc_10 sec instead of 13CPW90_10 μs). This test is only used on systems having a low band amplifier and probe capable of this rf field strength (typically 800 MHz systems).

¹³C RF homogeneity at limited power (15 microsecond)

Experiment - The power level is limited to an upper limit of the user-supplied value. A power sufficient to produce a ~15 s pw90 is used and the ¹³C pulse width is incremented. The data are plotted in a horizontal stacked array to allow inspection for arcing or phase instability. Based on the relevant maxima, a new array of 0, 360 and 720 degree ¹³C pulses is used to more accurately determine the rf homogeneity. The results are reported and stored in a separate history file from the standard history file (13CRFHOMOc).

¹³C RF homogeneity at limited power (10 microsecond)

Experiment - The power level is limited to an upper limit of the user-supplied value. A power sufficient to produce a ~10 μs pw90 is used and the ¹³C pulse width is incremented. The data are plotted in a horizontal stacked array to allow inspection for arcing or phase instability. Based on the relevant maxima, a new array of 0, 360 and 720 degree ¹³C pulses is used to more accurately determine the rf homogeneity. The results are reported and stored in a separate history file from the standard history file (13CRFHOMOc_10μs). This test is only used on systems having a low band amplifier and probe capable of this rf field strength (typically 800 MHz systems).

¹³C RF homogeneity at limited power (pw90=15.0 microsecond)

Experiment - The power level is limited to an upper limit of

the user-supplied value. The attenuator and modulator powers found for $\text{pw90}=15.0\ \mu\text{s}$ are used and the ^{13}C pulse width is incremented. The data are plotted in a horizontal stacked array to allow inspection for arcing or phase instability. Based on the relevant maxima, a new array of 0, 360, and 720 degree ^{13}C pulses is used to more accurately determine the rf homogeneity. The results are reported and stored in a separate history file from the standard history file (13CRFHOMOfc_15 μs).

^{13}C RF homogeneity at limited power ($\text{pw90}=10.0$ microsecond)

Experiment - The power level is limited to an upper limit of the user-supplied value. The attenuator and modulator powers found for $\text{pw90}=10.0\ \mu\text{s}$ are used and the ^{13}C pulse width is incremented. The data are plotted in a horizontal stacked array to allow inspection for arcing or phase instability. Based on the relevant maxima, a new array of 0, 360 and 720 degree ^{13}C pulses is used to more accurately determine the rf homogeneity. The results are reported and stored in a separate history file from the standard history file (13CRFHOMOfc_10 μs). This test is only used on systems having a low band amplifier and probe capable of this rf field strength (typically 800 MHz systems).

^{13}C RF amplifier linearity at limited power ($\text{pw90}=15$ microsecond)

Experiment - The ^{13}C pw360 is determined indirectly as a function of attenuator setting for ^{13}C on channel 2. The data are fitted to a linear regression. The maximum power is the attenuator setting found to produce a 15 μs pw90 .

Purpose - The amplitude response should be logarithmic. A log regression analysis should show the extent of fit to the ideal. The results are reported and stored in the history 13CAMPc.

^{13}C RF amplifier linearity at limited power ($\text{pw90}=10$ microsecond)

Experiment - The ^{13}C pw360 is determined indirectly as a function of attenuator setting for ^{13}C on channel 2. The data are fitted to a linear regression. The maximum power is the attenuator setting found to produce a 10 μs pw90 .

Purpose - The amplitude response should be logarithmic. A log regression analysis should show the extent of fit to the

ideal. The results are reported and stored in the history
13CAMPc_10 μ s.

$^{13}\text{C}/^{15}\text{N}$ power-limited decoupling tests

$^{13}\text{C}/^{15}\text{N}$ decoupling noise (FID) at user-selected decoupling power

Experiment - The time-domain noise characteristics in a single-scan experiment with no observe pulse is measured with:

- 1 no decoupling
- 2 ^{13}C decoupling
- 3 ^{15}N decoupling
- 4 $^{13}\text{C}/^{15}\text{N}$ decoupling, with power levels for decoupling set at the user-supplied limits.

The combined $^{13}\text{C}/^{15}\text{N}$ decoupling has 3 dB reduction for each channel from the user-supplied maxima. ^{15}N decoupling tests are skipped if the user-supplied maximum ^{15}N pulse power is set to zero. Average RMS noise, real/imaginary RMS, and real/imaginary dc offset results are stored in separate history files (NOISEc, NOISECc, NOISENc, and NOISECNc). Real and imaginary FIDs are plotted. Spectra are plotted.

Purpose - Decoupling should add minimal noise. This test permits a direct measurement of time-domain noise arising from decoupling at maximum permitted decoupling powers.

$^{13}\text{C}/^{15}\text{N}$ decoupling noise (FID) as function of decoupling power

Experiment - The time-domain average RMS noise in a single-scan experiment with no observe pulse is measured with:

- 1 no decoupling
- 2 ^{13}C decoupling
- 3 ^{15}N decoupling
- 4 $^{13}\text{C}/^{15}\text{N}$ decoupling, with power levels for decoupling varied from 1 dB up to the user-supplied limits.

The combined $^{13}\text{C}/^{15}\text{N}$ decoupling power has a 3 dB reduction for each channel from the user-supplied maxima. ^{15}N decoupling tests are skipped if the user-supplied maximum ^{15}N pulse power is set to zero. Average RMS noise is measured at each power level and the average of all these is stored along with the corresponding effective “loss” of signal-to-noise ratio (expressed in % loss). The average RMS

noise level and corresponding power levels are stored in a text file as the experiment proceeds. At the end, a histogram of the results is displayed and (optionally) plotted. The results are stored in separate history files (NOISE_c, NOISE_C_c, NOISE_N_c, and NOISE_CN_c)

Purpose - This test permits a direct measurement of the time-domain noise arising from decoupling as a function of power.

Sensitivity as function of ^{13}C decoupling power

Experiment - A single-scan experiment with a 90° degree observe pulse is measured with no decoupling and ^{13}C decoupling for decoupling power levels varied from 1 dB up to the user-supplied limit.

Signal-to-noise of the water is measured at each power level and the average of all these is stored along with the corresponding effective “loss” of signal-to-noise ratio relative to the no-decoupling results (expressed in % loss). The average signal-to-noise level and corresponding power levels are stored in a text file as the experiment proceeds. At the end, a histogram of the results is displayed and (optionally) plotted. A fixed region of noise is plotted for each power level to permit visual comparison. Plotted data also include printouts of signal-to-noise measurements and line widths for each spectrum.

The test is performed separately for cw, hardware modulator-based waltz-16, and waveform generator-based waltz-16 modulation. Results are stored in the SN_13Cdec_*_c history files.

Purpose - This test permits a direct measurement of frequency-domain noise arising from decoupling as a function of power, for different modulation schemes involving different hardware.

Sensitivity as function of ^{13}C decoupling power in presence of ^{15}N decoupling

Experiment - A single-scan experiment with a 90° observe pulse is measured with no decoupling; and ^{13}C decoupling with power levels for decoupling varied from 1 dB up to the user-supplied limit less 6 dB. ^{15}N decoupling at the user-supplied limit less 3 dB is applied in all experiments.

Signal-to-noise of the water is measured at each power level and the average of all these is stored along with the corresponding effective “loss” of signal-to-noise ratio relative to the no-decoupling results (expressed in % loss). The average signal-to-noise level and corresponding power levels are stored in a text file as the experiment proceeds. At the end, a histogram of the results is displayed and (optionally) plotted.

A fixed region of noise is plotted for each power level to permit visual comparison. Plotted data also have printouts of signal-to-noise measurements and line widths for each spectrum.

The test is performed separately for cw, hardware modulator-based waltz-16, and waveform generator-based waltz-16 modulation. Results are stored in the SN_¹⁵New_13Cdec_*_c history files.

Purpose - This test permits a direct measurement of noise arising from simultaneous ¹³C and ¹⁵N decoupling as a function of power, for different modulation schemes.

Sensitivity as function of ¹⁵N decoupling power

Experiment - A single-scan experiment with a 90° observe pulse is measured with no decoupling and ¹⁵N decoupling for decoupling power levels varied from 1 dB up to the user-supplied limit.

Signal-to-noise of the water is measured at each power level and the average of all these is stored along with the corresponding effective “loss” of signal-to-noise ratio relative to the no-decoupling results (expressed in % loss). The average signal-to-noise level and corresponding power levels are stored in a text file as the experiment proceeds. At the end, a histogram of the results is displayed and (optionally) plotted. A fixed region of noise is plotted for each power level to permit visual comparison. Plotted data also have printouts of signal-to-noise measurements and line widths for each spectrum.

The test is performed separately for cw, hardware modulator-based waltz-16, and waveform generator-based waltz-16 modulation. Results are stored in the SN_¹⁵Ndec_*_c history files.

Purpose - This test permits a direct measurement of noise arising from decoupling as a function of power, for different

modulation schemes involving different hardware.

Sensitivity as function of ^{15}N decoupling power in presence of ^{13}C decoupling

Experiment - A single-scan experiment with a 90° observe pulse is measured with no decoupling; and ^{15}N decoupling with power levels for decoupling varied from 1 dB up to the user-supplied limit less 6 dB. ^{13}C decoupling at the user-supplied limit less 3 dB is applied in all experiments.

Signal-to-noise of the water is measured at each power level and the average of all these is stored along with the corresponding effective “loss” of signal-to-noise ratio relative to the no-decoupling results (expressed in % loss). The average signal-to-noise level and corresponding power levels are stored in a text file as the experiment proceeds. At the end, a histogram of the results is displayed and (optionally) plotted.

A fixed region of noise is plotted for each power level to permit visual comparison. Plotted data also have printouts of signal-to-noise measurements and line widths for each spectrum.

The test is performed separately for cw, hardware modulator-based waltz-16, and waveform generator-based waltz-16 modulation. Results are stored in the SN_13Ccw_ $^{15}\text{Ndec}_*_\text{c}$ history files.

Purpose - This test permits a direct measurement of noise arising from simultaneous ^{13}C and ^{15}N decoupling as a function of power, for different modulation schemes.

^{13}C decoupling using phase modulation

^{13}C power level is reduced 20 dB from the hard-pulse level, limited by the user-supplied upper limit, and the ^{13}C decoupling efficiency is determined for the following phase-modulated, constant-amplitude broadband decoupling sequences using the built-in channel 2 hardware modulator (not the waveform generator).

The ^{13}C decoupling frequency is varied over a range of ± 80 ppm in a series of single-scan, proton-observe experiments. Only the ^{13}C -bound protons are shown in the expanded spectrum, which is plotted with spectra side-by-side in

absolute intensity mode to illustrate decoupling efficiency.

Purpose - This tests the hardware modulator under the user-supplied upper limit for decoupling.

^{13}C decoupling using adiabatic decoupling

The decoupling profile experiment is repeated using the (optional) channel 2 waveform generator under user-supplied upper limit on power. The waveforms are created automatically based on the calibrations produced in the AutoTest run and the user-supplied upper limit on power.

Purpose - This tests the waveform generation under the user-supplied upper limit for decoupling.

Decoupler heating

The same test as described above is performed, subject to the user-supplied upper limit on decoupling power. Results are stored in the history file `DECHEATc`.

Amplitude stability in the presence of ^{13}C decoupling

The 90° stability test is performed in the presence of ^{13}C decoupling at the user-supplied maximum decoupling power.

Purpose - Decoupling should not degrade amplitude stability which could lead to degradation of indirect detection experiments.

Amplitude stability in the presence of ^{15}N decoupling

The 90° stability test is performed in the presence of ^{15}N decoupling at the user-supplied maximum decoupling power. If the user-supplied maximum ^{15}N pulse power is set to zero, this test is skipped.

Purpose - Decoupling should not degrade amplitude stability which could lead to degradation of indirect detection experiments.

Amplitude stability in the presence of combined ^{15}N and ^{13}C decoupling

The 90° stability test is performed in the presence of the combined ^{15}N and ^{13}C decoupling at the user-supplied maximum decoupling powers less 3 dB each.

Purpose - Decoupling should not degrade amplitude stability which could lead to degradation of indirect detection experiments.

Decoupled methanol amplitude stability using a 6 kHz ¹³C decoupling RF field

The ¹³C-decoupled methanol signal is recorded under the user-supplied upper limit for ¹³C decoupling power. The ¹³C decoupling power is set to achieve a 6 kHz rf amplitude for 100 msec and a 1 second relaxation delay using waltz-16 modulation. The experiment is done both on- and off-resonance in ¹³C frequency. The off-resonance frequency is determined automatically to be at the edge of the waltz-16 effective bandwidth. Amplitude stability is measured for both and stored in a history file.

Purpose - It is important to characterize both on- and off-resonance decoupling since most experiments are designed for a wide variety of ¹³C chemical shifts. Incomplete or unstable decoupling would be most damaging at the extremes of decoupling bandwidths. This test measures this effect.

Decoupled methanol amplitude stability using a 3kHz ¹³C decoupling RF field

The ¹³C-decoupled methanol signal is recorded under the user-supplied upper limit for ¹³C decoupling power. The ¹³C decoupling power is set to achieve a 3 kHz rf amplitude for 100 msec and a 1 second relaxation delay using waltz-16 modulation. The experiment is done both on- and off-resonance in ¹³C frequency. The off-resonance frequency is determined automatically to be at the edge of the waltz-16 effective bandwidth. Amplitude stability is measured for both and stored in a history file.

Purpose - It is important to characterize both on- and off-resonance decoupling since most experiments are designed for a wide variety of ¹³C chemical shifts. Incomplete or unstable decoupling would be most damaging at the extremes of decoupling bandwidths. This test measures this effect.

Comparison of single-pulse versus HSQC decoupled methanol

Experiment - Single-scan experiments with a 90° observe

pulse are measured with ^{13}C decoupling. The same acquisition conditions are used except that a phase-cycled ^{13}C HSQC pulse sequence is used. The intensities and signal-to-noise ratios for each are averaged and the standard deviations calculated. The ratio of the average peak intensities and ratio of the average signal-to-noise ratios are calculated. The results are stored in the `13CSNc` history file.

Purpose - This test permits a direct measurement of the impact that the probe's rf homogeneity and the use of a multi-pulse sequence has in determining effective s/n in indirect detection experiments. In most probes the ^1H observe coil will detect signal from parts of the sample above and below the coil "window". This signal contributes to s/n in a single-pulse experiment, but does not in an X-edited experiment. The HSQC peak height is determined by, first, the portion of the sample accessible by the X-rf field (nominally the "window" size) and, second, by the X-coil rf homogeneity itself.

Comparison of single-pulse versus HSQC decoupled methanol (^{13}C pw90=10 microsecond)

Experiment - Single-scan experiments with a 90° observe pulse are measured with ^{13}C decoupling. The same acquisition conditions are used except that a phase-cycled ^{13}C HSQC pulse sequence is used. The intensities and signal-to-noise ratios for each are averaged and standard deviations calculated. The ratio of the average peak intensities and ratio of the average signal-to-noise ratios are calculated. The results are stored in the `13CSNc_10μs` history file.

Purpose - This test permits a direct measurement of the impact that the probe's rf homogeneity and the use of a multi-pulse sequence has in determining effective s/n in indirect detection experiments. In most probes the ^1H observe coil will detect signal from parts of the sample above and below the coil "window". This signal contributes to s/n in a single-pulse experiment, but does not in an X-edited experiment. The HSQC peak height is determined by, first, the portion of the sample accessible by the X-rf field (nominally the "window" size) and, second, by the X-coil rf homogeneity itself.

^{15}N Power-Limited Pulse and Decoupling Tests

If the standard AutoTest sample is augmented with $\sim 0.1\%$ acetonitrile (^{15}N , 99%), several tests analogous to those above are available for testing ^{15}N performance (using an indirect or triple-resonance probe). The tests are available for either channel 2 or channel 3 (with appropriate cabling). Results are stored in appropriate history files.

^{15}N pw90 and lowband amplifier compression

Experiment - The user-limited ^{15}N attenuator value is used. At this power level the pulse width is varied to obtain a pw90 value using an HMQC pulse sequence. Amplifier compression is determined by lowering power by 12 dB and re-determining the 90° pulse width.

^{15}N RF amplifier linearity at limited power

Experiment - The ^{15}N pw90 is determined indirectly as a function of the attenuator setting for ^{15}N .

The data are fitted to a linear regression.

Purpose - The amplitude response should be logarithmic. A log regression analysis should show the extent of fit to the ideal.

Acetonitrile ^{15}N -HMQC amplitude stability

Experiment - The ^{15}N HMQC acetonitrile signal reproducibility is recorded in multiple trials under the user-supplied upper limit for ^{15}N pulse power.

Purpose - ^{15}N indirect detection needs stable ^{15}N pulses to minimize ± 1 noise in 2D/3D experiments. This test measures the stability of the ^{15}N -edited spectrum.

Acetonitrile ^{15}N -HMQC decoupler heating

Experiment - The ^{15}N HMQC acetonitrile signal is recorded under the user-supplied upper limit for ^{15}N pulse power with 75 msec of decoupling at 2 KHz (scaled down for lower field systems) or the user-supplied upper limit for ^{15}N decoupling power. The temperature rise is measured. Plotted

data shows the speed to equilibrium and the temperature homogeneity during the process.

Purpose - ^{15}N indirect detection needs stable temperatures within the sample to minimize t_1 noise in multidimensional NMR experiments. This experiment measures the effect under typical conditions.

2H Pulse and Decoupling Tests

The standard AutoTest sample has 99% D_2O , which can be directly observed using the lock circuitry in the probe. Setting $\text{tn}=1\text{k}$ causes channel 2 to produce ^2H rf, which is then routed into the ^2H Diplexer box and then on to the lock BNC of the probe. This setup is used for ^2H gradient shimming. Since ^2H is observed, the pw90 can be determined at any power level. The normal channel 2 hardware setup is used, and because ^2H is a low-frequency nucleus, the rf is automatically routed through the first broadband amplifier before going to the rf relay at the magnet leg. If tn is not 1k , the rf is routed to the broadband preamp for normal observe. With suitable control of rf from a decoupling channel, the ^2H pulses can be delivered from a decoupling channel, and thus calibrated. The pw90 s determined can be used to calculate dmf/dpwr values for heteronuclear decoupling experiments on the channel used. Results are stored in appropriate History files. Usually, only one particular channel is used for ^2H decoupling, so only one of the tests below is used for a particular spectrometer.

^2H pw90 using channel 2

Experiment - The power level of channel one's attenuator is varied for a fixed pw of $150\ \mu\text{s}$. At the power giving a maximum signal, the pulse width is then varied to obtain a calibrated pw90 . The power level and pw90 values are stored in the history file H2PW90 .

^2H pw90 using channel 3

Experiment - The power level of channel three's attenuator is varied for a fixed pw of $150\ \mu\text{s}$. At the power giving a maximum signal, the pulse width is varied to obtain a calibrated pw90 . The power level and pw90 values are stored in the history file H2PW90_ch3 .

^2H pw90 using channel 4

Experiment - The power level of channel four's attenuator is varied for a fixed pw of 150 μs . At the power giving a maximum signal, the pulse width is varied to obtain a calibrated pw90. The power level and pw90 values are stored in the history file H2PW90_ch4.

^2H pw90 using lock/decoupler

Experiment - The power level of the Lock/Decoupler's attenuator is varied for a fixed pw of 150 μs . At the power giving a maximum signal, the pulse width is varied to obtain a calibrated pw90. The power level and pw90 values are stored in the history file H2PW90_1kdec. The Lock/Decoupler, if present, is normally the last channel.

^2H 90 degree pulse amplitude stability and sensitivity

Experiment - The ^2H D₂O signal is recorded following 90° pulses on channel 1. Average amplitude and standard deviation are recorded and stored in the H2STAB90 history file. Sensitivity and linewidth are measured and stored in the H2SENSITIVITY history file.

Purpose - ^2H heteronuclear decoupling involves 90° pulses and constant amplitude, phase-modulated pulse trains. A stable amplitude is essential for optimum results. High sensitivity gives a more stable lock. This test permits tracking of lock sensitivity over time.

^2H 90 degree pulse amplitude stability using channel 3

Experiment - The ^2H D₂O signal is recorded following 90° pulses on channel 3. Average amplitude and standard deviation are recorded and stored in the H2STAB90ch3 history file.

Purpose - ^2H heteronuclear decoupling involves 90° pulses and constant amplitude, phase-modulated pulse trains. Stable amplitude is essential for optimum results.

^2H 90 degree pulse amplitude stability using channel 4

Experiment - The ^2H D₂O signal is recorded following 90° pulses on channel 4. Average amplitude and standard

deviation are recorded and stored in the H2STAB90ch4 history file.

Purpose - ^2H heteronuclear decoupling involves 90° pulses and constant amplitude, phase-modulated pulse trains. Stable amplitude is essential for optimum results.

^2H 90 degree pulse amplitude stability using lock/decoupler

Experiment - The ^2H D_2O signal is recorded following 90° pulses using a Lock/Decoupler. Average amplitude and standard deviation are recorded and stored in the H2STAB90lkdec history file.

Purpose - ^2H heteronuclear decoupling involves 90° pulses and constant amplitude, phase-modulated pulse trains. Stable amplitude is essential for optimum results.

^2H spinlock test

Experiment - The ^2H D_2O signal is recorded following a spinlock pulse on channel 1. Average amplitude and standard deviation are recorded and stored in the H2SPNLKSTAB history file.

Purpose - ^2H heteronuclear decoupling involves 90° pulses and constant amplitude, phase-modulated pulse trains. The lock circuitry should be able to provide a high-power spinlock.

Installation Tests for Cryogenic Probes

After a cryogenic probe is installed a series of experiments is performed to both condition the probe and to verify its operation. These experiments run a ^{13}C HSQC pulse sequence with combinations of ^{13}C and ^{15}N decoupling during acquisition. The experiments run for 15, 30, 60, or 90 minutes followed by a test in which the power of ^{13}C and ^{15}N decoupling is varied while noise in the FID is measured (no observe pulse is used). The results are plotted to determine the extent of any excess noise from decoupling. The test is usually repeated automatically until stopped by the user.

This test conditions the probe by heating the coil enough to drive off (mainly) absorbed water. This water vapor (in the vacuum chamber of the probe) can give plasma noise under an rf field which adds to the noise of the FID and degrades s/n ratios. Per the probe conditions, the residual water is pumped away by the vacuum pump.

15-minute conditioning (ATcryo15)

Experiment - 15 minutes of conditioning is done followed by the noise test.

30-minute conditioning (ATcryo30)

Experiment - 30 minutes of conditioning is done followed by the noise test.

60-minute conditioning (ATcryo60)

Experiment - 60 minutes of conditioning is done followed by the noise test.

90-minute conditioning (ATcryo90)

Experiment - 90 minutes of conditioning is done followed by the noise test.

Other Test Descriptions

This section describes additional tests in AutoTest.

Heating during spin lock test

Experiment - The same test as in the variable temperature test is performed using a 70 ms ^1H pulse at an rf field strength of 10 kHz with a total recycle time of 1.5 s, including acquisition. The data are plotted in the same way as in the VT test.

Purpose - Many modern experiments use spin locks or decoupling within their pulse sequences. Often, rf fields can cause significant sample heating. Depending on the nature of the probe, this heating can be a problem causing baseline artifacts and t_1 noise. It is important to quantify the amount of sample heating, the speed in attaining a new equilibrium temperature, the homogeneity of temperature throughout the sample during the heating period, and the final change in temperature.

This test imposes a rather strong (10 kHz) rf field for a period of time often used in TOCSY experiments, using a recycle time of 1.5 s (including acquisition). Single transients

are acquired, at a rate of one per 1.5 s. The data show any temperature change at the expected 0.01 ppm/degree. The intensity and/or linewidth can be used to measure temperature homogeneity. The number of transients needed to attain a new equilibrium temperature measures the ability of the probe to stabilize the effects of internal sample heating. The final shift value indicates the total temperature change. This can be used to reduce the requested temperature value so as to obtain the desired equilibrium. Of course, the amount of heating is a function of the sample itself, primarily its salt content.

The same approach may be used to follow the actual temperature in the sample under the influence of X-nucleus decoupling.

Lock tests

Experiment - Lock power is varied over a 30 dB range and the lock level recorded. The experiment is repeated for the lock gain. A log regression analysis is performed to confirm the relationship between the lock signal and power/gain.

Purpose of Lock Gain Test - Lock gain is selectable in a logarithmic manner (in dB). In an ideal case, variation of receiver gain should produce a logarithmic dependence of signal strength.

Purpose of Lock Power Control Using an RF Attenuator - Overall lock power control is accomplished using computer-controlled rf attenuators. The amplitude response should also be logarithmic. A log regression analysis should show the extent of fit to the ideal.

Spectral purity test

Experiment - Four single-scan 100 kHz spectral width spectra are acquired with no pulse. Each spectrum is plotted with a few millimeters of noise.

Purpose - RF purity of the transmitter and receiver can be tested by recording data without any excitation pulse. The spectrum reveals any artifactual signals.

Variable temperature test

Experiment - Single-scan spectra are acquired during an increase of 5° C in sample temperature. Spectra are recorded sequentially. Spectra are taken every 2 s until the

sample reaches equilibrium, as reflected in a stable chemical shift of a methyl proton. Spectra are plotted in a stacked manner to permit examination of the rate of change of temperature, the homogeneity of temperature, and the length of time necessary to reach equilibrium.

Purpose - The sample temperature is increased by 5° C while recording spectra every 2 s. Most methyl resonances show a chemical shift of $(\text{sfrq}/100)$ Hz/° C and this shift, therefore, indicates the actual temperature distribution within the sample. The methyl resonance should move quickly and homogeneously to its new equilibrium position. The rate of change and homogeneity of change demonstrate the VT performance of the probe and regulation hardware.

Small-angle phase shift test

Experiment - Single-scan spectra are acquired in which the phase of the pulse is incremented in each spectrum through a full 360° at constant receiver phase. Spectra are plotted in a horizontal stack to show a smooth phase rotation of the spectrum.

Purpose - Small-angle phase adjustment is used in multiple-quantum selection ($q > 2$), phase-modulated pulses, and a variety of complex pulse sequences. This test exercises the phase-shift hardware by varying the pulse phase in small increments over 360°.

Tests Using Salty Sample

The following tests are designed for a standard AutoTest sample with added salt.

¹H RF homogeneity

Experiment - One hundred experiments are run in which the pulse width is incremented from 1 to 100 s. The spectra are plotted in a horizontal stack in phased mode and sufficiently expanded so that the base of the water can be examined using the same phase settings for each spectrum.

Pulse stability

Experiment - A single-scan pulse experiment is repeated 20 times and the spectra plotted in a horizontal stack. The average peak amplitude, rms deviation and sensitivity are

measured and reported. This test is run for the following:

- 90° flip pulses
- 30° flip pulses
- 10° flip pulses
- 1 μ s pulses

Purpose - The same tests are performed as for the standard AutoTest sample, but for a sample with added salt. Any difference in performance is highly relevant for user samples containing salts or buffers, as these constituents detune the probe and result in lower sensitivity and longer pw90's.

Heating during spin lock test

Experiment - This test is performed using a 70 ms ^1H pulse at an rf field strength of 10 kHz with a total recycle time of 1.5 s, including acquisition. The presence of salt provides opportunity for additional sample heating and this test measures the extent of heating under the same conditions as used for the non-salty sample.

Decoupler heating

The same test as in the variable temperature test is performed but this time using a 75 ms ^{13}C decoupling period prior to acquisition within a total recycle time of 1.5 s, including acquisition. One-hundred, single-scan spectra are collected with ^{13}C decoupling followed by 100 identical spectra with no decoupling. The spectra are plotted in a stacked manner to permit examination of the rate of change of temperature, the homogeneity of temperature, and the length of time necessary to reach equilibrium. The rf field strength is sufficient to decouple over a 160 ppm range using garp-1. The presence of salt provides opportunity for additional sample heating and this test measures the extent of heating under the same conditions as used for the non-salty sample.



Agilent Technologies

© Agilent Technologies, Inc.

Printed in USA, March 2013



G7446-33301